

# DermaPlast® SOFT

## General Product Description/Intended Purpose

DermaPlast® SOFT strips and woundplasters are non-active, single-use, self-adhesive, non-sterile devices consisting of an adhesive coated backing material (carrier material) and a woundpad, protected with a silicone paper.

The main constituents of this products are:

- The backing / carrier material:
  - skin color non-woven
- The medical grade adhesive:
  - Transparent synthetic rubber based adhesive / Hotmelt
- The wound pad with 2 different layers:
  - Exudate-transmission layer made of viscose.
  - Absorption layer of non-woven with absorbent fibers.
- Protective paper / film for the protection of the wound pad and adhesive side of the product:
  - Silicone paper

The products do not contain resorbable components, medicinal products / drug substances (derived from human blood or plasma or not), tissues or cells of human or animal origin, radioactive substances, nanomaterials or other hazardous components.

The strips are sealed individually in peel-pouches and packed in cardboard folding boxes. Wound plasters are packed directly in cardboard folding boxes.

## Application/Indication

DermaPlast® SOFT strips and woundplasters are non-active, single-use, self-adhesive, non-sterile devices consisting of an adhesive coated backing material (carrier material) and a woundpad, protected with a silicone paper, for the protection of minor injuries and absorption of exudate.

The products are intended for humans, without restrictions (in age group, weight range, health status or condition), for use by lay-users and healthcare professionals, in contact with intact skin (healthy skin) and breached skin, and short-term (60 minutes ≤ Duration ≤ 30 days, acc. to MDR) / prolonged use (24 hours < Duration ≤ 30 days, acc. to ISO 10993-1:2018) duration of use.

## Reference Numbers

Packaging	Dimensions	Name	Reference number
P1	4cmx5m	DermaPlast SOFT	535351
P1	6cmx5m	DermaPlast SOFT	535361
P1	8cmx5m	DermaPlast SOFT	535371
P200	22mm	DermaPlast SOFT Spots	535382

# DermaPlast® SOFT

## Contra-Indications

Not indicated to be applied on postsurgical wounds, highly bleeding wounds, ulcers or chronic wounds.

## Single Use Device

Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

## Product Disposal

To minimize the risk of potential infection hazards, or environmental pollution, disposable components should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

## Incident Reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

## Product Performance Characteristics

Property	Test method	Unit	Specifications
Visual examination	EN ISO 5350-02	Conform/not conform	Conform
Peel test (90° peel angle)	EN ISO 7500-01:2016-05	N/25mm	17 +3 -5
Absorption capacity	(DIN) EN 13726-1:2002:06	g/100cm <sup>2</sup>	>6

## Labelling

### Lot-No. with 8-Digit Code

e.g.:



0

12

XXXXX

year

week of production

for internal purposes only

### Manufacturing Date

e.g.:



2020

03

01

year

month

day

### Expiry Date

e.g.:



2025

03

01

year

month

day

Shelf Life: 5 years

# DermaPlast® SOFT

Medical Device



Distributor Symbol



**Latest Date of Revision:** 2020-03-26