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As a specialist for fitness and health products, Ludwig Artzt GmbH has been bridging the gap between sports science and the world of fitness for more than 35 years. Resistance devices constitute a particular focal point of our product range. Under the ARTZT vitality® umbrella brand, specialists and home users have access to a comprehensive selection of high-quality products and modern, efficient training concepts for professionals and recreational athletes.

Ongoing, close communication with experts is the pillar that supports the success of Ludwig Artzt GmbH. Our own training and continuing education institute is considered one of the key points of contact for qualified seminars on the topics of health and fitness. Disseminating technical knowledge backed up by sports science is the focus of our educational activities.

The elastic bands can be found exclusively in many exercise programs. Clinical research has proven that elastic resistive exercise increases/improves:

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Material and characteristics

- » pure latex* with low protein content
- » superior resistance
- » linear build-up of force that is easy on the joints and prevents injuries
- » easy to transport and extremely versatile in application
- » available in five colour-coded levels of resistance (yellow / red / green / blue / black)
- » each in three lengths (2.5 m / 6.0 m / 25.0 m) for individual training options

Widerstand/resistance	
leicht/easy	
mittel/medium	
stark/strong	
extra stark/extra strong	
spezial stark/specially strong	

*

  - Medical product contains pure natural latex, may cause allergic reactions.
- Avoid enduring contact with skin and mucosal areas.

Factors that influence durability and safe exercising with the bands:

- » Protect the bands by keeping it away from sharp objects. Finger rings, long, sharp fingernails and rough soles of sport shoes could damage the bands. Always examine the bands before use for small nicks, tears, or punctures that may cause the band to break. The use of the additional accessories will make exercising easier and protects the bands.
- » Store the bands out of direct sunlight and away from extreme temperatures.
- » It is best to keep it rolled up or hanging.
- » If the bands become sticky, clean with mild soap and water and dry flat.
- » Use of the bands can cause serious injuries when not used properly or safe. Therefore always perform safe and painfree exercises. Avoid exercises that involve stretching the band in such a fashion that it may snap toward the head and cause injuries to the head or eyes. If these types of exercises are prescribed, protective eyewear should be worn.
- » The exercise bands are not toys. Supervision is recommended when children are exercising with the bands.
- » Avoid choosing bands that are too short. Bands should be between 2.5 m and 3.0 m long to avoid dangerous exercises and damage of the band.

Don't overstretch the band by more than 3 times its resting length (for example, 30 cm to more than 90 cm).

- » After use in chlorinated water, rinse the bands with tap water and dry flat.

Selection of the correct band


The appropriate colour / resistance of the band depends on your individual physical condition. Choose a band that allows you to complete 15 repetitions.

Experience collected during the use of elastic bands in practice has shown that yellow and red exercise bands are mainly used in rehabilitation after injuries and for seniors. The green, blue and black exercise bands from ARTZT vitality® are especially suitable for use in sports and fitness.

Wrapping and connecting of the bands

Your elastic band should be securely attached to your hand or foot before use to avoid slippage and possible injuries. Double-wrapping the band may help to secure it to your hand or foot.

Attachelement to the hands


Wrapping at one hand

Lay the band flat in your hand with the end towards your pinky finger. Wrap the long end of the band around the back of your hand. Repeat as needed. Grasp.

Wrapping with both hands


Begin with palms up and end of band between the thumb and finger. Rotate your palms inward, bringing the band around the back of your hands. Repeat as needed. Grasp.

Attachment to the feet

Foot Loop

Loop over the top of the foot and stabilize other ends with opposite foot.

Foot Wrap

Stand in the middle of the band. Wrap one end around the top of the foot.

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The main aims of a health-oriented preventive exercise program are increasing endurance, muscle mass and to improve bodyshape and posture. To match those goals, beginners should choose 6 to 8 exercises from the following program. Complete 1 to 3 sets of 15 to 20 repetitions with mild fatigue on the last set. At no time should you feel the strain of the exercise as very heavy; if so, change to an easier band colour. You may feel more comfortable using different levels of resistance bands for different body parts! Progress to the next colour of band or add sets or training units to your program when you feel comfortable and you are able to easily complete the 3 sets of 15 to 20 repetitions. Follow these guidelines when exercising with bands:

- » Perform motions slowly and under control. Maintain good overall body posture before, during and after each movement.
- » The band should be under slight tension before each movement to ensure the correct joint position.
- » The direction of the pull of the band influences both muscles and joints. Therefore it is important to be aware of the „force angle“ (angle between the band and the lever arm (arm, leg)). The maximum torque

(strength) is often noted at 90° (mid-range). Below 30° there is little significant effect.

- » Do not exercise while experiencing pain. If you feel pain stop exercising. If your pain should persist consult a health care provider.
- » All exercises are shown in the starting and final position. Closely examine the photos and imagine the movement. Always exercise BOTH sides of the body!

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Stand on the middle of the band with one foot slightly in front of other.

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Elbows bent and close to the body. Palms up. a) Bend elbows and slowly return.

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b) Variation „backhoe„ Same exercise as a) but move your arms up while bending your elbows. (go on with lunge)

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Stand on the middle of the band with one foot slightly in front of other. Elbows are close to the body and the forearms are rotated outward.

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a) Lift arms to the side up to shoulder level (sliding movement!) and slowly return.

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b) Stand in parallel position on the middle of the band. Cross ends of bands in front of your body. Lift arms to the side up to shoulder level and slowly return.

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Stand in parallel position. Begin with band wrapped around upper back, holding each end in hands. Elbows bent.

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Push arms forward, straightening elbows. Hands are lightly rotated inward and finally form clenched fists.

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Stand in parallel position. Grasp the mid of the band (distance between the hands about 10 cm) at shoulder level with arms extended.

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Keep one arm extended and pull the other arm back. Change arms. Shoulder stays low, the elbow high. Relax arms after 2 to 3 repetitions and start again.

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Stand on the middle of the band with one foot and grasp the ends of the band at the level of your navel. The elbows are close to the body and the palms face each other.

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a) Extend both arms upward and to the outside

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b) Variation Airman: Step with your back foot to the front. Lean forward and extend both arms and the back leg against band.

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Stand in the middle of the band. Cross end of band in front of your body. Keep palms up and elbows to side. The right elbow is bent with exterior rotation while the left hand is resting loosely on the hip.

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Now extend the right leg and the left arm simultaneously. Switch sides.

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Lie on back with hips and knees flexed. Stretch middle of band around both feet and wrap ends of band around each foot. Grasp both ends of band and slightly put the hands behind the neck.

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a) Curl trunk upwards. The head remains in the extension of the torso while the gaze is directed towards the ceiling at an angle. Hold briefly and then lower your upper body. Ideal for beginners, because the band supports the exercise.

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b) With the help of a partner, the band comes from behind. Slightly raise your upper body and pull your extended arms in the direction of your thighs. The band creates resistance during the exercise.

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Lie on back with the middle of the band under your shoulder blades. The legs are bent and the soles of the feet touch the floor. The right arm is bent and the upper arms point up. The left arm is resting loosely on the floor with the band wrapped around the hand.

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Alternating between right and left arm, lift upper body and pull band along the opposite thigh.

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V-seat with the legs slightly bent. Wrap the band around the sole of the right foot and grasp the ends of the band with both hands (fists) in front of the chest.

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Now slowly rotate the upper body to the left against the resistance of the band. Hold briefly and then change sides (feet).

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Wrap middle of band around sole of one foot and grasp the ends of the band with both hands.

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Keeping back and neck straight, extend leg backwards against band. Slowly retrace. A soft pad (for example Stability Trainer) will minimize the pressure on the floor knee.

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Wrap the two ends of the band around the left foot once and stand on the ends of the band. The right foot is inside the "small" loop of the band.

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Extend the right leg by lifting the right foot up to the side. Change sides.

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Wrap the two ends of the band around the left foot once and stand on the ends of the band. The right foot – offset slightly towards the front – is inside the "small" loop of the band.

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Now move the right leg to the left, passing in front of the supporting leg. Lead the movement with the heel! Change sides.

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Stand with one foot slightly in front of other. Wrap middle of band around the back foot and step with the front foot on both bands. Grasp the ends of the band in front of the body with both hands.

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Flex knee, lifting band upward. Slowly return.

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Step on the middle of the band and grasp the ends of the band in front of the hip area. Now bend the knees and move your buttocks BACK and down as though you wanted to sit down on a chair. Raise the arms to shoulder height at the same time.

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Extend the legs while simultaneously stretching both arms upwards. The gaze follows the hands during this process.

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Grasp the middle of the band in front of your body. One hand at hip height the other about 20 cm beside. The left hand holds the band in front of the body with the elbow bent.

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a) Bend down while reaching to the front and down with the left hand (picking up a ball). Then return to the initial position.

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b) Extend the left arm up and to the left after standing up from the bent position. Follow the movement with your eyes. Combine a) + b) several times. Switch sides.

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Stand in parallel position. Wrap the band around your thighs and secure the ends of the band in place.

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a) Now take several long side-steps to the left. Change sides. You can also combine steps to the right and left.

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b) During the side-step, bend down and pick up an imaginary object

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Stand on the middle of the band. Grasp ends of band with tension beside your thighs.

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a) Bend your knees and slightly come back to the straight position against the resistance of the band. Please note, that the band stays slightly tensed during the exercise.

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b) Put the band on your shoulders and bend and extend your knees.

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The family builds the basis for a healthy physical development of children. Therefore it is very important that parents exercise with their kids in leisure time. This could be easily done at home if it is raining outside. Just try it!

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Wrap the band around the waist and start. Pull your partner through the room or just move against the resistance of the band. Is it possible to explore the room with eyes closed?

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Mum and dad stretch the band around 10 cm above the floor. Balance on the band (forward, backward, with eyes closed...).

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Stand on the middle of the band. Grasp both ends of the band in front of your body. The band is slightly tensed.
a) Now the fight starts. Box in the air, against an invisible opponent, against dad's stomach.

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b) Variation Superman / Superwoman – Grasp both ends of the band beside your body and lift your hands to the ceiling (single-arm or bilateral).

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b) Put different coloured bands on the floor and hop over the „zebra crossing“.

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

Mum and dad stretch bands in different heights and diagonals. How could I get on the other side of the labyrinth? (crawling, jumping...)

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©2018 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.
Ludwig Artzt GmbH - Schiesheck 5 - 65599 Dornburg - Germany
©2018 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.