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## Usage and Care Instructions Dynair® Pro

**Congratulations on your purchase of the Dynair Pro.**

**Having purchased the Dynair Pro you have chosen a TOGU branded product which is “Made in Germany”!**

**It is made in conformity with the applicable EU regulations for sports equipment, it contains no heavy metals and is made of 100% recyclable plastic.**

The first ball cushion combined with a stable base - two pieces of equipment rolled into one. Providing even more training options for sensorimotor functions and coordination. Particularly suitable for balance exercises on one leg - can also be used upside down! With extra skid-proof material for absolute safety. The Dynair Pro can be used in pairs or in rows depending on the sport (e.g. football, golf, tennis, skiing, etc.).

### Warnings:

- The screws in the base plate are correctly tightened at the factory. Please do not retighten them!
- Do not store the Dynair Pro in direct sunlight or near sources of heat as these may cause the product to go out of shape.
- Avoid contact with sharp and pointed objects (e.g. glass splinters / stones in soles of shoes) as these can damage the Dynair Pro. Do not throw the product!
- Wear sports shoes / golf shoes (but not with metal studs) when using the Dynair Pro or train barefoot but never use the product while wearing socks.
- Check the Dynair Pro for signs of damage before use as the product cannot be guaranteed safe for use if damaged.
- The Dynair Pro must be placed on a horizontal, stable and dry floor surface to ensure that it cannot skid off to one side during use.

### Apparatus data:

The Dynair Pro is made from high-quality Rulon and has an inflation valve.

Dimensions: approx. Ø disc 36 cm x 18 cm height

Load capacity: approx. 200 kg

### Care:

Any mild standard cleaning agents and disinfectants can be used to clean the Dynair Pro. Solvent-based, acidic or caustic cleaners will damage the material. The Dynair Pro is waterproof. First and foremost, you should always make sure that the soft anti-skid coating (circular rubber surface on base of Dynair Pro), which gives the Dynair Pro its exceptional stability, is always clean and dry.

### Inflation:

If necessary, the air pressure in the Dynair Pro can be adjusted with the pump supplied. The Dynair Pro is fitted with a needle valve for this purpose. Do not change the internal pressure of the Dynair Pro (inflation or deflation) until the Dynair Pro has reached room temperature and do not inflate the Dynair Pro to the extent that the surface is no longer flat. Otherwise the product may be overstretched or insufficiently tensioned.

Always moisten the pin on the pump before inserting it in the needle valve.

### Storage:

Please do not store your Dynair Pro upright.