promed EMT-4

DDX model: DX6605BO

Digitales, elektrisches Schmerztherapie- und Muskelstimulations-Gerät (TENS & EMS)
(TENS: Transkutane Elektrische Nerven-Stimulation) (EMS: Elektrische Muskel-Stimulation)

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




PRODUCT INFORMATION / ACCESSORIES

About this product:

Your **Promed EMT-4** is a battery-operated unit with two adjustable output channels. This combined TENS-EMS unit generates electrical pulses whose amplitude, duration and modulation can be varied via the buttons. The **Promed EMT-4** can be prescribed by your doctor (classification number 09.31.02.1-999). The Promed EMT-4 device is only intended for domestic use.

- 1** Output sockets for electrode cables
- 2** On/Off button
- 3** Program select button
- 4** Intensity regulation buttons
- 5** Confirm button
- 6** Channel select button
- 7** Battery compartment with 4 x 1,5 V batteries, AAA type
- 8** Electrode plugs (cable – electrodes)
- 9** Connector plugs (cable – main unit)

Product features:

- 1** Eleven preset programs (see program table) and two programmable programs
- 2** LCD (Liquid Crystal Display) for: **a)** Treatment time; **b)** Program; **c)** Frequency (pulses per second); **d)** Pulse width (duration of individual pulses); **e)** Intensity
- 3** Two channels, independently adjustable in 50 stages
- 4** Connection check: when the electrodes are not connected or attached, the intensity cannot be set higher than stage 1, or reverts to "0". (to avoid draining the batteries)
- 5** Application protection: when the unit is switched on, the intensity is set to stage 0 (of 50).
- 6** Memory function: the unit stores the treatment duration.
- 7** Time display: displays the time in "sleep mode".

PACKAGE CONTENTS

System components

Your **Promed EMT-4** unit includes the following accessories::

- 1 TENS/EMS main unit: 1
- 2 Electrode connector cables: 2
- 3 Self-adhesive electrodes: 2 pairs (4 pads)
- 4 Unit stand: 1
- 5 Operating instructions: 1
- 6 Batteries: 4
- 7 Storage case: 1

INTRODUCTION

What is pain?

Pain is the body's own early warning system. Feeling pain is important, as it indicates an abnormal condition within the body and serves as a warning to us before any further damage or injury occurs. However, long-term, persistent pains – often known as chronic pain – serve no obvious useful purpose once a diagnosis has been made. TENS was developed to soothe or eliminate certain types of chronic and acute pain.

D

US

F

I

E

INTRODUCTION

We differentiate between two fundamental types of pain:

- **Acute pain**

as a chief symptom can often help the doctor with diagnosis and the acute pain has a protective function for the patient.

- **Chronic pain**

can often become part of the illness in itself. A patient suffering from chronic pain will often suffer for years and experiences changes in his/her personality structure.

The term “stimulation current therapy” serves as an umbrella term for the electrical stimulation of nerves and muscles.

Basic terms in electrical pain therapy

Pulse frequency

Number of pulses per second, measured in Hertz (Hz).

Tens works between 1 and 200 Hz.

Pulse width

The duration of an individual pulse, measured in microseconds.

1µsec. corresponds to 1/1.000.000 second.

Current intensity

Measured in milliamperes 1mA = 1/1000 A.

Tens units work with a maximum of 50 mA.

Modulation

Pulse intensity, pulse frequency and pulse duration are simultaneously varied.

INTRODUCTION

How does TENS work?

TENS is a non-invasive and drug-free method of treating pain. The TENS unit transmits soothing electrical pulses through the skin that stimulate the nerve (or nerves) in the area being treated. In most cases, this stimulation greatly reduces or even completely eliminates the sensation of pain, as it blocks the original pain signals to the brain.

It is also reported that TENS stimulation promotes the release of endorphins into the bloodstream, leading to a further reduction in pain.

The benefits of TENS devices in pain management have been clinically proven in the case of many patients. To get the most out of your TENS unit, read these instructions carefully and follow your doctor's recommendations for treatment.

All information, general sensations, skin sensations, deep sensibility, is relayed to the CNS (central nervous system) via the nerve pathways.

The classic TENS method **interrupts pain transmission** for an indefinite period of time.

It does not affect the source of pain locally, but acts **exclusively on the nerve tissue**.

Nerve tissue can be stimulated by electrical current **over a wide spectrum of frequencies**.

Because of the interaction between the stimulation frequency and the transmission speed of the nerve fibers, it can be seen **that different frequencies selectively activate certain types of nerve**.

D

US

F

I

E

INTRODUCTION

The body has two different types of nerve fiber:

- a) Unmyelinated nerve fibers with continuous, slow relaying of stimuli.
- b) Myelinated nerve fibers with saltatory, fast relaying of stimuli.

MOTOR (myelinated) NERVE FIBERS (skeletal muscles)

Up to 10 Hz	Single twitches
10 – 20 Hz	Fibrillations
20 – 80 Hz	Tetanic contractions
above 90 Hz	Muscle waves

AUTONOMIC (unmyelinated)

NERVE FIBERS (smooth musculature, organs)

Below 10 Hz	Stimulation of the vascular muscles
10 – 20 Hz	Stimulation of the intestinal muscles*

* smooth musculature (muscles of the internal organs, deep lying – involuntarily contracting muscles)

INTRODUCTION

Treatment theories

1. Endorphin theory

At very low frequencies ($< 10\text{Hz}$), stimulation current predominantly causes hormonal manipulation of vegetative centers in the central nervous system. Pain-transmitting fibers are stimulated, triggering central pain-inhibiting mechanisms and opioid release in the brain stem or diencephalon, and the body's own endorphins and morphine-like substances are released. These work like painkillers!

2. Gate control theory

Stimulation frequencies $> 10\text{Hz}$ have their effect mainly in the "dorsal horns" of the spinal cord. There, pain sensations are prevented from penetrating to consciousness. Frequencies from 85 – 100 Hz to 200 – 250 Hz can trigger pain-inhibiting mechanisms in the free nerve endings and sympathetic nerve fibers.

This is the most widespread theory.

General pain analysis

1. Intense, sharp or piercing pain

Use frequencies $< 50\text{ Hz}$ to $< 20\text{ Hz}$

Select as high a pulse width as possible $> 150\ \mu\text{sec}$

2. Dull, throbbing pain

Use frequencies $> 50\text{ Hz}$ (typically)

Select as wide a pulse width as possible (the pain or throbbing must not be aggravated), that usually means a small pulse width of $< 120\ \mu\text{s}$

Avoid stimulating muscles at the same time, otherwise this can cause aching muscles!

D


US

F

I

E

INTRODUCTION

How TENS works

Tens electrical pulses override the pulses from the source of pain.

POSSIBLE TENS APPLICATIONS

Scientifically founded, TENS proves effective in treating the following pain conditions:

1. Acute pain conditions:

- Postoperative pain
- Fracture pains
- Childbirth pains

2. Chronic pain conditions:

- Pain in the neck, back and small of the back
- Painful stiffness in the shoulders
- Rheumatic complaints
- Ischialgia
- Arthrosis joint pains
- Tendovaginitis
- Pain from Bekhterev's arthritis
- Extremity pains

* PC1 or PC2:

120 Hz, select pulse width (μS)

according to the part of the body:

- Face $50 \mu\text{S}$
- Shoulder/upper arm $150 \mu\text{S}$
- Lower arm $200 \mu\text{S}$
- Thighs $350 \mu\text{S}$
- Lower leg/foot $400 \mu\text{S}$

D

US

F

I

E

POSSIBLE TENS APPLICATIONS

Generally the Promed EMT-4 shall be used for the following medical indications or the following complaints as a TENS unit for treatment:

- For symptomatic relief of chronic intractable pain
- For post traumatic pain (acute pain occurring)
- For post surgical pain (pain caused by surgery)

It is imperative to consult a doctor BEFOREHAND in the following cases:

- Prolonged pain symptoms despite therapy
- When using strong painkillers or local anesthesia
- Infectious diseases
- Circulation problems, thrombosis
- Impaired sensibility (numbness)
- Use with babies and infants
- Pregnancy
- Psychosis
- Bleeding tendency
- Cancer
- Extreme sensitivity to or fear of electricity
- Pacemaker
- Artificial joints or metal plates

TENS THERAPIE RECOMMENDATIONS

Therapy goal	Frequency	Alternative
Skin stimulation	approx. 50 Hz	approx. 1–4 Hz
Analgesia (pain relief)	approx. 50–180 Hz	approx. 2–4 Hz
Detoning (relaxing, loosening)	approx. 100–120 Hz	approx. 2–8 Hz
Trophic (tissue nutrition)	approx. 100–120 Hz	

Indications:

- for symptomatic relief of chronic intractable pain
- for post traumatic pain (acute pain occurring)
- for post surgical pain (pain caused by surgery)

HOW DOES EMS (ELECTRICAL MUSCLE STIMULATION) WORK?

EMS stands for **E**lectronic **M**uscle **S**timulation.

Muscle contractions are triggered by electrodes positioned on the body. The successive contracting and relaxing of the muscles strengthens them and at the same time helps integrate the muscles more quickly into general movement patterns, promoting and stimulating both muscle activity and arterial circulation.

The **Promed EMT-4** is ideally suited to preventing and delaying disuse atrophy in injured parts of the body, e.g. when a part of the body is immobilized for a long period of time due to a plaster cast, splint or extension bandage.

A further area of application for the **Promed EMT-4** is the treatment of various mobility disorders, for example involving the hip, ankle and elbow joints. The stimulator can be used on the corresponding parts of the body to improve coordination in the case of mobility disorders. This acts on both the **nerve** and

muscle tissue. By loosening the muscles, pain can be relieved or even eliminated. **The stimulation frequencies are lower than with the TENS method.**

During the entire cycle time, the muscles being treated can be

- **tensed**
- **relaxed**
- **moved**
- **kneaded**
- **massaged**

thereby achieving a clear reduction in local pain.

Both the TENS and EMS methods can activate the muscles and thus have a muscle-building effect.

D

US

F

I

E

HOW DOES EMS (ELECTRICAL MUSCLE STIMULATION) WORK?

Low frequencies (1-10Hz)

with long pulse times and single muscle contractions provide:

- detoxification
- relaxation
- improved circulation in the muscles
- carrying away metabolic end-products (lymph drainage)

Medium frequencies (20-50Hz)

The quick consecutive contractions (fibrillation) provide:

- increased strain on the muscles
- increased muscle building
- increased muscle fiber size
- improved oxygen uptake
- increased energy consumption for the user

POSSIBLE EMS APPLICATIONS

Generally the Promed EMT-4 shall be used for the following medical indications or the following complaints as a EMS unit for treatment:

Indications:

- For relaxation of muscle spasm
- For increase of blood circulation
- For prevention of retardation of disuse atrophy
- For muscle re-education

POSSIBLE EMS APPLICATIONS

Applications

- Regenerating parts of the body that have been at rest for long periods (plaster cast)
- Treating other mobility disorders, e.g. with the hip, ankle and elbow joints
- Supporting training performance and muscle building
- Warming up the muscles
- Muscle stamina training
- Muscle relaxation
- Muscle detox
- Improving circulation and metabolism
- Rehabilitation after injuries
- Pain therapy
- Stimulating arterial circulation
- Promoting muscle activity, especially in older people
- Improving coordination in the case of mobility disorders
- Preventing thrombosis
- Treating scoliosis
- Cellulitis/toning: this shows how versatile **Promed** electrotherapy can be. With regular use, combined with an exercise and nutritional program, you can achieve visible success. All TENS/EMS pulse current types are suitable here. Several treatments a day are recommended.

D

US

F

I

E

MASTER POINTS AS A SOURCE FOR IMPROVING PERFORMANCE


These points, stemming from acupuncture, can be used at any time as support or can even be chosen as main points. The master points are to be found on the lower arm, chest and leg.

Lower arm master points (Pe 6, between the tendons of the radial flexor muscles of the wrist on the inside of the lower arm, 2 Cun proximal to the wrist and Lu 9, on the radial side of the fold between hand and wrist, lateral to the Arteria Radialis) have a relaxing effect on the cardiovascular system, and improve circulation and blood flow properties.

Chest master points (Lu 1, one Cun below the collar bone, six Cun lateral to the center line between the first and second ribs) increase breathing capacity.

Leg master points (Ma 36, one finger's breadth lateral to the lower edge of the tuberositas tibia/tibial plateau, three Cun below the knee fold) are said to be the best general tonifying points.

The master points can be used individually or in combination. They are present on both the left and right sides of the body.

CONTRAINDICATIONS

If you apply to one of the following cases always be sure to consult a doctor before you apply the Promed EMT-4 and clarify the use of the device with him.

- In case of persistent pain symptoms, despite therapy
- Taking a strong analgesic or local anesthetic
- In Infectious Diseases
- By circulatory disorder (thrombosis and embolism)
- For sensory disturbances (numbness)
- When used in infants, toddlers and children
- In pregnancy
- With Psychosis
- If bleeding tendency
- With Cancer
- At extreme sensitivity to or fear of electricity
- In patients with metallic implants
- For heart problems, especially cardiac arrhythmias
- Prior to each electrode placement, that apply current to the carotid sinus region (anterior neck)
- Prior to each electrode placement, that apply current transcebrally (through the head)
- If undiagnosed pain symptoms
- For treatment of the eyelids
- Serious arterial circulatory problems (embolism) in the lower limbs

- The device should not be used for symptomatic local pain relief unless etiology is established or unless a pain syndrome has been diagnosed
- When cancerous lesions are present in the treatment area
- At swollen, infected, inflamed areas or skin eruption (e.g. phlebitis, thrombophlebitis, varicose veins, etc.)
- On-existence of a demand-type cardiac pacemaker or any implanted defibrillator
- In areas of the body with poorly enervated areas
- With epilepsy
- For abdominal or inguinal hernia

WARNINGS

- Pregnant women should not be treated with Promed EMT-4 during the first trimester.
- TENS devices must be kept out of the reach of children.
- The safety of using TENS devices during pregnancy or delivery has not been proven.
- TENS is not effective in treating pain of the central nervous system (headaches).
- Should TENS treatment fail to work or cause discomfort, stimulation. Should be discontinued until a doctor has reassessed the situation.

D

US

F

I

E

WARNINGS

- Always switch the TENS unit off before attaching or removing the electrodes.
- Never position the electrodes over the eyes, in the mouth or within the body.
- TENS devices do not possess any healing properties.
- TENS is a method for treating symptoms and as such only suppresses the perception of pains that would otherwise serve as a defense and warning mechanism for the body.
- Do not apply stimulation while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk of injury.
- If you are in the care of a physician, consult with your physician before using this device.
- If you have had medical or physical treatment for your pain, consult with your physician before using this device.
- If your pain does not improve, becomes more than mild, or continues for more than five days, stop using the device and consult with your physician.
- Do not apply stimulation over your neck because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Do not apply stimulation across your chest because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.
- Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Do not apply stimulation when in the bath or shower.
- Do not apply stimulation while sleeping.
- Do not use the device on children, if it has not been evaluated for pediatric use.
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Apply stimulation only to normal, intact, clean, healthy skin.
- In the application you need to be careful when in the immediate vicinity a mobile phone is switched on, as this may affect the output of the device.
- If you are pregnant, you should not place any electrodes directly over the uterus or on the belly. The reason is that the heart of the foetus may be affected (so far there are no medical reports of damage).
- The user may be the intended operator.
- Do not modify the device or the electrodes without authorization of the manufacturer. This could cause improper functioning.
- Because the effects of stimulation on the brain are unknown, do not use the device on opposite sides of your head.
- If you had surgery recently, consult your physician before you use the device.
- Using the device may disrupt the healing process.
- This device is designed for use by and on a single adult person. For hygienic reasons electrodes should not be shared.

WARNINGS

- Do not modify electrodes (e.g. by cutting them). This increases the current density, which is potentially hazardous an effective current density beyond 2 mA/cm² requires increased awareness).

PRECAUTIONARY MEASURES

- TENS is not a substitute for pain medications and other pain management therapies.
- TENS devices have no curative value.
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- The success of the therapy depends to a considerable extent on the patient's choice of physician and whether they are qualified for the treatment of patients suffering from pain.
- The long-term effects of electrical stimulation are unknown.
- You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium.
- Use caution if you have a tendency to bleed internally, such as following an injury or fracture.

PRECAUTIONARY MEASURES

- Consult with your physician prior to using the device after a recent surgical procedure, because stimulation may disrupt the healing process.
- Use caution if stimulation is applied over the menstruating or pregnant uterus.
- Use caution if stimulation is applied over areas of skin that lack normal sensation.
- Use this device only with the leads, electrodes, and accessories recommended by the manufacturer.
- Make sure you do not use the unit if it or the cables are damaged or faulty. If a fault occurs, please send in the unit or cables via your specialist dealer to:
Promed GmbH
Service Department,
Lindenweg 11,
D-82490 Farchant.

D


US





F

I

E


POSITIONING ELECTRODES

- D
- US**
- F
- I
- E

DISPLAY GUIDE

When the unit is switched on, all readouts are active for a short time before the unit switches to standby (wait) mode.

PROGRAM SELECTION...

Preset programs P-01 to P-11 and memory function

1. After attaching the electrodes, press the „ON/OFF” button to switch on the unit.
2. Keep pressing the “Program” button until one of the preset programs (P-01 to P-11) is shown on the display.
3. Press the “Set” button. The time will flash on the display, “ ’ ” for minutes and “ ` ” for seconds. Use the “+” or “-” buttons to enter the desired treatment duration (20 minutes to 2 hours) in steps of ten minutes.
4. Press the “Select” button 2 times. The intensity display “CH1” will start to flash. Use the “+” or “-” buttons to enter the desired value (from 1 to 50). Setting the intensity will simultaneously start the unit. In the top right of the display, the clock symbol will begin to flash. Press the “Select” button once more. The intensity display “CH2” will start to flash. Once again, use the “+” or “-” buttons to set the desired value.

AND STARTING THE TREATMENT

5. To complete the settings, press the "Select" button once more.
6. If you want to change the intensity during treatment, press the "Select" button again and adjust the intensity as required using the "+" or "-" buttons.

MEMORY FUNCTION

Calling up and resetting the memory

1. Press the "ON/OFF" button to switch on the unit.
2. Press and hold down the "Set" button for more than 3 secs. to switch to the memory function. The total treatment time will now be displayed for 30 secs.
After this, the unit will switch back to standby mode, or you can press any button to return to setting mode.
3. When in memory mode, hold down the "Select" button for 3 secs. to reset the treatment time to zero.

If no buttons are pressed for at least 2 minutes before confirming a function, the unit will automatically return to time display mode. To switch off the unit, press the "ON/OFF" button once more.

D

US

F

I

E

PROGRAM SELECTION AND STARTING THE TREATMENT

Programmable programs P-C1 and P-C2

After attaching the electrodes...

1. Press the "ON/OFF" button to switch on the unit.
2. Keep pressing the "Program" button until the programmable programs (P-C1 or P-C2) are shown on the display.
3. Press the "Set" button once. In the display, the frequency (Hz) and "SET" will start to flash. Use the "+" or "-" buttons to enter the desired value (1 Hz to 200 Hz).
4. Press the "Set" button. In the display, the pulse width (μS) and "SET" will start to flash. Use the "+" or "-" buttons to enter the desired value (from 30 μS to 400 μS) in steps of ten.
5. Press the "Set" button once more. The time will flash on the display. Use the "+" or "-" buttons to enter the desired treatment duration (1 minute to 24 hours).
6. Press the "Set" button to confirm and then press the "Select" button 2 times. The intensity display "CH1" will start to flash. Use the "+" or "-" buttons to enter the desired value (from 1 to 50). Setting the intensity will simultaneously start the unit. Press the "Select" button once more. The intensity display "CH2" will start to flash. Once again, use the "+" or "-" buttons to set the desired value. Fast adjustment protection will be activated and stops at intensity level 15 when the intensity increases rapidly.
7. To complete the settings, press the "Select" button once more.
8. If you want to change the intensity during treatment, press the "Select" button again and adjust the intensity as required using the "+" or "-" buttons.

PROGRAM SELECTION AND STARTING THE TREATMENT

Pause mode

1. To temporarily pause treatment, press the "Program" button. In the display, the pause symbol will start to flash. To continue, press the "Program" button once more.
2. To stop treatment ahead of time, press the "ON/OFF" button 2 times. The unit will revert to standby mode and display the time.
3. Once the treatment period is completed (countdown), the unit will automatically switch to standby mode.

Setting the time

1. In standby mode, the unit display the time. If it is showing the treatment settings, press the "ON/OFF" button.
2. To set the clock, now press the "Set" button. The hours will start to flash. Use the "+" or "-" buttons to set the desired hour. Press the "Set" button once more. The minutes will start to flash. Use the "+" or "-" buttons to set the minutes. Press the "Program" button to confirm.

D

US

F

I

E

PROMED EMT-4 PROGRAM TABLE

Programm	Pulse 1		Pulse 2		TENS/EMS	Quantity per cycle	Which program for which purpose?
	Frequency (Hz)	Pulse width (µS)	Frequency (Hz)	Pulse width (µS)			
P01	15	200	Burst frequency = 2 Hz		TENS	Burst, 54 times, 56 secs. break 5 secs. per cycle	Extremities (legs, arms) and back, warming up, loosening up, relaxing
P02	150	30	Burst frequency = 2 Hz		TENS	Burst, 54 times, 29 secs. break 5 secs. per cycle	Face, headaches, promotes concentration
P03	5	100	Burst frequency = 2 Hz		EMS	Burst, 3 times, 5,8 secs. break 11,6 secs. per cycle	Loosening up, stamina, aerobic metabolism
P04	75	100	60	300	EMS	Modulation, 5 secs. break 10 secs. per cycle	Muscle building, hypertrophy, anaerobic metabolism
P05	50	175	-	-	EMS	Conventional, 5 secs. break 10 secs. per cycle	Strength endurance, anaerobic metabolism, rehabilitation
P06	110	175	-	-	EMS	Conventional, 5 secs. break 10 secs. per cycle	Arms, maximum strength, intramuscular coordination
P07	100	125	100	225	EMS	Modulation, 5 secs. break 10 secs. per cycle	Arms, legs, maximum strength, intramuscular coordination
P08	130	100	-	-	TENS	Conventional, 31 times break 5 secs. per cycle	Extremities, back, shoulders, pain relief
P09	35	30	Burst frequency = 2 Hz		TENS	Burst, 40 times, 29 secs. break 5 secs. per cycle	Neck, face, shoulders, pain relief, loosening up
P10	40	100	-	-	TENS	Conventional, 62 secs. break 5 secs. per cycle	Shoulder pain, knee complaints
P11	50	250	-	-	EMS	Conventional, 5 secs. break 10 secs. per cycle	Breast toning, pectoralis training, rehabilitation
PC1	1-200	30-400	-	-	TENS	-	Individual treatment, free to be set and stored by patient or doctor
PC2	1-200	30-400	-	-	TENS	-	Individual treatment, free to be set and stored by patient or doctor

TECHNICAL SPECIFICATIONS

- Power supply: 4 x batteries 1,5 V, AAA type
(total 6,0 V) or 4 x rechargeable batteries 1,2 V AAA type
- 2 channels: CH 1 and CH 2
- Current type: biphasic square pulse
- Intensity: 0 mA – 50 mA, at 500
(each channel is separately adjustable, 50 stages)
- Frequency range: 1 Hz to 200 Hz
- Pulse width range: 30 μ sec to 400 μ sec
- treatment time:
 - a) Preset programs: (P-01 to P-11):
20 mins to 2 hours
 - b) Programmable programs: (P- C1 and P- C2):
1 min to 24 hours
- Environmental conditions:
 - Temperature: 5 °C to 40 °C
 - Relative humidity: 30 to 75 %
- Storage conditions:
 - Temperature: - 10 °C to + 55 °C
 - Relative humidity: 10 to 90 %
- Dimensions: 122 x 58 x 18 mm (L x W x H)
- Service life of the device: 3 years
- Service life of battery: New batteries will last for approx. one week (when used for once a day, P08, 50 mA intensity)
- Service life of electrode pads: Pads can be cleaned and reused for up to 10 -15 times
- Applied part: Electrode
- IP classification system: IP22

D

US

F

I


E

FURTHER INFORMATION

Certification No.

Classification label

V
mA
mW
Hz
 μ s
kOhm
mm
g
BF
P

Volts
Milliamps
Milliwatts
Hertz
Micro second
Kilo-Ohms
Millimetres
Grammes
Body floating
Output rating

General warning.
According to DIN VDE 750/219, the unit must be marked with this symbol. It means that the user must follow these operating instructions and thus use the unit for the purpose intended.

Direct Current

TROUBLESHOOTING

If the Promed EMT-4 is not functioning correctly

1. Make sure that the batteries are correctly inserted, or replace the batteries. When replacing the batteries, make sure that the polarity is correct.
2. If the display is showing all functions but you can not feel any stimulation, check that the electrode cables are properly connected and the electrodes are correctly positioned. If the unit seems to be working, but you feel no stimulation, it is possible that the electrode cables or the electrodes have to be replaced.
3. If the unit seems to be working, but you still feel no stimulation, increase the current intensity as described above.
4. If you still can't feel any stimulation, the electrodes or electrode cables may have to be replaced.

SERVICE AND MAINTENANCE

Should you encounter further problems, please contact your dealer or return the unit to your dealer. Do not attempt to repair a faulty unit yourself.

For repair or maintenance send the **Promed EMT-4** to the following address via your specialist dealer:

Promed GmbH

Service Department

Bischofstraße 1, D-82490 Farchant

Tel. +49 (0) 8821 96 21-0, Fax +49 (0) 8821 96 21-21

D

US

F

I

E

DISPOSAL

Power tools, accessories, and packaging should be recycled in an environmentally suitable manner.

Only for EU countries:

Do not throw power tools in with the domestic waste! In accordance with the European Directive 2012/19/EU on waste electrical and electronic equipment and its implementation into national law, nonfunctioning power tools must be collected separately and recycled in an environmentally suitable manner.

This symbol indicates that disposal of this product in the domestic waste is not permitted within the EU. Waste devices contain valuable recyclable material that should be recycled. Recycling also prevents uncontrolled waste disposal from damaging the environment and human health. Therefore, please dispose of all waste devices using the appropriate collection systems or send the device to the place of purchase for disposal. The place of purchase will then recycle the device.

Rechargeable batteries/batteries:

Do not throw your rechargeable batteries/batteries in the domestic waste, in fire, or in water. Rechargeable batteries/batteries should be collected, recycled, or disposed of in an environmentally friendly manner.

Only for EU countries:

According to the directive 91/157/EEC, defect or used rechargeable batteries/batteries must be recycled. Waste rechargeable batteries/batteries can be directly handed to:

Promed GmbH, Lindenweg 11, D-82490 Farchant.

GUARANTEE CONDITIONS AND CUSTOMER SERVICE

This product was produced with the utmost care and was thoroughly tested before it left the factory. Therefore, we provide a warranty period of 24 months from the date of purchase for this product. **Promed** products conform to their description and the relevant specifications; it is your responsibility to ensure that the product which you have purchased is designed for your intended use.

- In the event of material defects and/or manufacturing errors which appear during proper use and are detected during the warranty period, we will replace all defective parts of the product free of charge within the warranty period, including the labour costs for warranty repairs.
- **The warranty does not cover:**
 - Normal wear on the product
 - Defects which were caused during transportation or storage of the product
 - Defects or damage which were caused through improper use or inadequate maintenance
 - Damage resulting from disregarding the notes in the operating manual
 - Damage resulting from modifications to the product which were not carried out by **Promed**.
 - Damage caused by sharp objects or as a result of torsion, compression, being dropped, an abnormal impact or other actions which are beyond **Promed's** reasonable control.
 - Parts which are subject to wear (e.g. moving parts such as ball bearings, etc., fasteners, etc.) are generally excluded from the warranty.
- **A claim under the warranty is not enforceable if:**
 - the product is not sent back in its original packaging or in adequately secure packaging;

GUARANTEE CONDITIONS AND CUSTOMER SERVICE

- it has been modified or repaired by a person or company other than **Promed** or one of **Promed's** authorised dealers;
- the product has been repaired using parts which are not approved by **Promed**;
- the serial number / LOT number of the product has been removed, deleted, changed or made illegible.
- For reasons of hygiene, any products which have been exposed to direct contact with a body or with bodily fluids (e.g. blood) must be packaged in an additional plastic bag before they are sent back. In this case, a specific note of this situation must be made in the package or in the accompanying documents.
- For components or products which are replaced during the warranty period, a warranty is only provided for the remainder of the original warranty period; provided that the product or part replacement was carried out by **Promed** or one of **Promed's** authorised dealers.
- The warranty period begins on the day of purchase. Claims under the warranty must be established within the warranty period. Complaints arising after the expiry of the warranty period cannot be considered.
 - Under these warranty conditions, the warranty shall only apply if the date of purchase can be proven by means of a sales receipt or similar.
 - We reserve the right to make technical and visual changes as well as changes to the features.
- This warranty is only valid and enforceable in the country in which it was purchased by the original buyer, provided that it was **Promed's** intention that the product should be offered for sale in this country. This warranty is also enforceable in any country within the European Economic Area in which **Promed** has an authorised importer or sales partner. Depending on each country, specific and differing warranties and guarantees may be applicable owing to the particular applicable legislation. These legal provisions are neither excluded nor limited by these warranty conditions. To the extent permissible under national law, the warranty period shall not be extended, renewed or otherwise affected by any subsequent resale, repair or replacement of the product.
- The provisions of the UN Convention on Contracts for the International Sale of Goods shall not apply.
- Our warranty conditions do not affect the seller's legal warranty obligation.
- To the fullest extent possible under the applicable mandatory legislation, these warranty conditions constitute your sole and exclusive legal remedy and shall apply in place of any other express or implied warranty conditions. **Promed** is not liable for any unusual, incidental, punitive or consequential damages including, but not limited to loss of profits, loss of use, loss of revenue, costs of replacement equipment or facilities, third party insurance claims, and damage to property which are the result of the purchase or use of the product or occur owing to a breach of warranty, a breach of contract, negligence, product defects or any other legal or statutory circumstances, even though **Promed** was aware of the possibility of such damages. **Promed** is not liable for any delay in warranty claims.
- **Promed** cannot be held liable for any possible translation errors.
- The following information is vital for problem-free processing:
 1. Original sales slip/receipt or dealer's stamp with the date of purchase
 2. Defect detected
 3. Unit name / type / serial / batch number

D

US

F

I

E

CLINICAL RECOMMENDATIONS

Clinical recommendations for Promed EMT-4 settings

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




E

ELECTROMAGNETIC COMPATIBILITY (EMC)

- This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used.
- Use of accessories other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Table 1

Declaration - electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for domestic establishment and in establishment directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2

Declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°		
	0 % UT; 250/300 cycles		
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

ELECTROMAGNETIC COMPATIBILITY (EMC)

Table 3

Declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3V 0.15 MHz to 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz	Not applicable	<p>Portable and mobile RF communications equipment should be used no closer to any part of device, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d=2.1 \sqrt{P}$ 150 KHz to 80 MHz $d=2.1 \sqrt{P}$ 80 MHz to 800 MHz $d=2.3 \sqrt{P}$ 80 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

- Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which device is used exceeds the applicable RF compliance level above, device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating device.
- Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and device, as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	0.15 MHz to 80 MHz $d=2.1 \sqrt{P}$	80 MHz to 800 MHz $d=2.1 \sqrt{P}$	80 MHz to 2.7 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer			
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Symbol for **WARNING!** Indicates sections of the user manual that are safety-relevant!

Instructions for a unit, system or accessories contain important safety-relevant information for the consumer.

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Symbol for Type **BF (body floating)** units

Appliances that comply with the requirements of this norm offer a higher level of protection against electric shocks, especially when taking into account the permissible patient leakage currents and safety of the user.

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Please note the symbol for instruction manual:

The instruction manual must be followed effectively.

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Symbol for **Manufacturer and manufacturing Factory**

It shows the name and the address of the manufacturer and the manufacturing factory.

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Symbol for date of manufacture
It shows the date of manufactur.

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Symbol for Notified Body
This Symbol Complies with the European Medical Device Directive 93/42/EC and amended by directive 2007/47/EC requirements. Notified body TÜV Rhineland (0197).

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Symbol for Authorized representative in the European Community
It shows the Name of Authorized representative in the European Community.

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Symbol for Authorized representative in the switzerland.

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To indicate the entity importing the medical device into the locale

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The first number 2: Protected against solid foreign objects of 12,5 mm and greater.

The second number: Protected against vertically falling water drops when enclosure tilted up to 15°. Vertically falling drops shall have no harmful effects when the enclosure is tilted at any angle up to 15°, on either side of the vertical.

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Shenzhen Dongdixin Technology Co., Ltd.
Shanghai International Holding Corp. GmbH (Europe)
MedNet SWISS GmbH, Bäderstrasse 18, 5400 Baden,