promed IT-6

MODEL: LT3006

Digitales, elektrisches Inkontinenztherapiegerät (TENS & EMS)
(TENS: Transkutane Elektrische Nerven-Stimulation) (EMS: Elektrische Muskel-Stimulation)

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GETTING STARTED	p. 25	NORMALIZED SYMBOLS	S. 239
CONTENTS	p. 26	CE DIRECTIVES	S. 243
SYSTEM COMPONENTS	p. 27	TECHNICAL SPECIFICATIONS	S. 245
DESCRIPTION OF THE UNIT	p. 28	TRANSPORT/STORAGE/OPERATING CONDITIONS	S. 247
INTRODUCTION	p. 29-31	ACCESSORIES	S. 269
INDICATIONS / CONTRAINDICATIONS	p. 31/32	CLINICAL RECOMANDATIONS	S. 273
WARNINGS	p. 32/33	GUARANTEE CARD	S. 278
CAUTIONS / CONTRAINDICATIONS	p. 33/34	MANUFACTURER + FURTHER INFORMATION	S. 0
OPERATING THE UNIT	p. 34-40		
PROGRAM	p. 41/42		
SAFETY INSTRUCTIONS	p. 43		
TROUBLESHOOTING	p. 44		
CARE, MAINTENANCE AND STORAGE	p. 45		
DISPOSAL	p. 45		
GUARANTEE CONDITIONS	p. 46		

DE

EN

FR

IT

ES

NL

RU


PL

SE

FI

Your **Promed IT-6 - package** contains the following equipment:

- | | |
|-----------------------------|-------|
| 1. IT-6 Unit | 1 PCS |
| 2. Cable | 2 PCS |
| 3. Electrode pads (40*40mm) | 4 PCS |
| 4. Vaginal probes | 1 PCS |
| 5. User manual | 1 PCS |
| 6. AAA battery | 4 PCS |
| 7. Carrying case | 1 PCS |
| 8. Belt Case | 1 PCS |

We reserve the right to make changes for the purpose of technical improvement without prior notice.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Please read the user manual thoroughly before using the device for the first time.

Thank you for purchasing the **Promed IT-6**. You have purchased a high-quality product designed for personal hygiene and wellbeing purposes. **Promed** is a leading company with decades of experience in personal hygiene, wellness and health areas. The **Promed IT-6** unit has been designed and manufactured in accordance with the Medical Device Directives 93/42/EEC and 2007/47/ EC to guarantee quality during use. The device may be used once this user manual has been read.

We, as the manufacturer, cannot be made liable in any way for injury or damages to people or objects that arise from failure to comply with this user manual. We wish you a lot of enjoyment with your new **Promed IT-6**. We would like to familiarize you with the therapie unit in the following sections. Please read the user manual thoroughly before using the device for the first time.

You can receive professional advice wherever **Promed** products are sold or get in touch with us if you have any questions. We can provide you with the name of the representative responsible for you.

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Further informations can be found on our homepage:

www.promed.de

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Promed IT-6 electrical stimulator is a portable electrotherapy device featuring two therapeutic modes: Transcutaneous Electrical Nerve Stimulation (TENS), and Electrical Muscle Stimulation (EMS), which are used for pain relief and electrical muscle stimulation. The device is used for the rehabilitation of incontinence, weakness of the perineal muscles and for pain relief. By using this device regularly, you can strengthen your pelvic floor muscles significantly, and the improvement of your muscle function will help you to live your life again free of stress, urge mixed or faecal incontinence. The stimulation device transmits via electrodes that are applied to the skin or lighten the electrical current to underlying nerves and muscles on the introduced vaginal-anal probe. The parameters of the device are controlled by the press buttons. The intensity level is adjustable according to the needs of patients. Before using, please read all the instructions in this user manual carefully and keep it safe for future reference.

What is pain?

Pain is the body's own early warning system. Feeling pain is important, as it indicates an abnormal condition within the body and serves as a warning to us before any further damage or injury occurs. However, long-term, persistent pains – often known as chronic pain – serve no obvious useful purpose once a diagnosis has been made. TENS was developed to soothe or eliminate certain types of chronic and acute pain.

We differentiate between two types of pain:

- **Acute pain**
as a chief symptom can often help the doctor with diagnosis and the acute pain has a protective function for the patient.
- **Chronic pain**
can often become part of the illness in itself. A patient suffering from chronic pain will often suffer for years and experiences changes in his/her personality structure.

What is TENS?

Transcutaneous Electrical Nerve Stimulation (TENS) is a non-invasive, drug free method of controlling pain. TENS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. TENS does not cure any physiological problem; it only helps control the pain. TENS does not work for everyone; however, in most patients it is effective in reducing or eliminating the pain, allowing for a return to normal activity.

How does TENS work?

Das **Promed IT-6** transmits harmless electrical signals into the body by the self-adhesive electrodes or the vaginal / anal probe. This relieves pain in two ways:

- Firstly, it blocks the body's pain signals which are normally transmitted from the area of damage through the nerve fibres to the brain - TENS interrupts these pain signals.
- Secondly, TENS stimulates the body's production of endorphins – its own natural painkillers.

Possible TENS applications

Generally the **Promed IT-6** shall be used for the following medical indications or the following complaints as a TENS unit for treatment:

- For symptomatic relief of chronic intractable pain
- For post traumatic pain (acute pain occurring)
- For post surgical pain (pain caused by surgery)

What is EMS?

Electrical Muscle Stimulation (EMS) is an internationally accepted and proven way of treating muscular injuries. It works by sending electronic pulses to the muscle needing treatment; this causes the muscle to exercise passively. It is a product derived from the square waveform, originally invented by John Faraday in 1831. Through the square wave pattern it is able to work directly on muscle motor neurons. This device has low frequency and in conjunction

with the square wave pattern allows direct work on muscle groupings. This is being widely used in hospitals and sports clinics for the treatment of muscular injuries and for the re-education of paralyzed muscles, to prevent atrophy in affected muscles and improving muscle tone and blood circulation.

How does TENS work?

The EMS units send comfortable impulses through the skin that stimulate the nerves in the treatment area. When the muscle receives this signal it contracts as if the brain has sent the signal itself. As the signal strength increases, the muscle flexes as in physical exercise. Then when the pulse ceases, the muscle relaxes and the cycle is repeated. The goal of electrical muscle stimulation is to achieve contractions or vibrations in the muscles. Normal muscular activity is controlled by the central and peripheral nervous systems, which transmit electrical signals to the muscles. EMS works similarly but uses an external source (the stimulator) with electrodes attached to the skin for transmitting electrical impulses into the body. The impulses stimulate the nerves to send signals to a specifically targeted muscle, which reacts by contracting, just as it does with normal muscular activity.

Possible EMS application

Generally the **Promed IT-6** shall be used for the following medical indications or the following complaints as a EMS unit for treatment:

- For relaxation of muscle spasm
- For increase of blood circulation
- For prevention of retardation of disuse atrophy
- For muscle re-education

What is INCONTINENCE?

Urinary incontinence, the involuntary loss of urine from the bladder, is a problem for many people. There are many types of incontinence. Urinary incontinence, which is divided into stress incontinence and urge incontinence, and the faecal incontinence, in which the involuntary sphincter cannot hold my bowels. The electrical stimulation through a vaginal / anal probe, or surface electrodes in some cases is a well-tolerated treatment for urge, stress, and mixed incontinence and has shown the positive results in improving bladder and bowel control of the patients.

Stress incontinence

Stress incontinence is urine loss by increased abdominal pressure on the bladder which is caused during coughing, sneezing, laughing, playing sports or lifting something heavy. Stress incontinence is the most common form of incontinence which primarily affects the women. Usually they occur if the perineal and pelvic muscles are weakened, such as through pregnancy, childbirth and menopause.

Urge incontinence

Urge incontinence is a sudden, strong urge which is followed by an immediate contraction of the bladder leading to an involuntary loss of urine. Both men and women may be affected thereof, especially this form of incontinence applies to the older people. A reason for this may be a fault in the part of a nervous system, which controls the function of the bladder.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Mixed incontinence

Mixed incontinence is a combination of stress incontinence and urge incontinence.

Anal incontinence

Anal incontinence also referred to as faecal incontinence is the impairment of the ability to control the outflow of gas or stool. There are many possible causes of anal incontinence, the most common is the injury of the anal sphincter (circular muscle), for example, caused during birth, by surgery or by damage to the nerves which controls the anal muscles. The condition usually worsens with age.

How does an incontinence treatment functions?

The electrical stimulation over the pelvic nerve is an acceptable alternative for the treatment of urinary incontinence. It is also recommended as a treatment method for anal incontinence which is caused by pelvic floor dysfunction or a malfunctioning of anal muscles. In the treatment of stress incontinence, the aim of the electrical stimulation is to reflect a voluntary muscle contraction, and thus to improve the function of pelvic floor muscles. For urge incontinence, the aim is to inhibit involuntary contractions of the bladder by stimulation of the nerves in the pelvic floor. In the treatment of mixed incontinence, stimulation is suitable to be used for both urge incontinence and stress incontinence. For anal incontinence the aim is to improve the bowel control by strengthening and toning pelvic floor muscles.

Be sure to consult a doctor before using the unit!

In the USA, the sale of this device is restricted by law to doctors, or subject to a doctor's prescription.

Please carefully read and understand the following warnings and cautions to ensure the safe and correct use of this device and to prevent injury.

Indications

TENS: The device is designed to be used for temporary relief of pain, including the acute and chronic pain relief. EMS: The device is designed to be used for stimulate muscles for prevention of muscular atrophy, increase of blood flow circulation, muscle strengthening and improve and facilitate muscle performance. For rehabilitation of incontinence, weakness of the perineal muscles.

**Contraindications**

Do not use this device without consulting a doctor if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death. If you apply to one of the following cases always be sure to consult a doctor before you apply the **Promed IT-6** and clarify the use of the device with him.

- In case of persistent pain symptoms, despite therapy
- Taking a strong analgesic or local anesthetic
- In Infectious Diseases
- By circulatory disorder (thrombosis and embolism)
- For sensory disturbances (numbness)

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

INDICATIONS / CONTRAINDICATIONS

- When used in infants, toddlers and children
- In pregnancy
- With Psychosis
- If bleeding tendency
- With Cancer
- At extreme sensitivity to or fear of electricity
- In patients with metallic implants
- For heart problems, especially cardiac arrhythmias
- Prior to each electrode placement, that apply current to the carotid sinus region (anterior neck)
- Prior to each electrode placement, that apply current transelebrally (through the head)
- If undiagnosed pain symptoms
- For treatment of the eyelids
- Serious arterial circulatory problems (embolism) in the lower limbs
- The device should not be used for symptomatic local pain relief unless etiology is established or unless a pain syndrome has been diagnosed
- When cancerous lesions are present in the treatment area
- At swollen, infected, inflamed areas or skin eruption (e.g. phlebitis, thrombophlebitis, varicose veins, etc.)
- On-existence of a demand-type cardiac pacemaker or any implanted defibrillator
- In areas of the body with poorly enervated areas
- With epilepsy
- For abdominal or inguinal hernia
- Do not apply stimulation over your neck, (especially the carotid sinus), because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure;

WARNINGS

General

- Pregnant women should not be treated with **Promed IT-6** during the first trimester.
- TENS devices must be kept out of the reach of children.
- The safety of using TENS devices during pregnancy or delivery has not been proven.
- TENS is not effective in treating pain of the central nervous system (headaches).
- Should TENS treatment fail to work or cause discomfort, stimulation should be discontinued until a doctor has reassessed the situation.
- Always switch the TENS unit off before attaching or removing the electrodes.
- Never position the electrodes over the eyes, in the mouth or within the body.
- TENS devices do not possess any healing properties.
- TENS is a method for treating symptoms and as such only suppresses the perception of pains that would otherwise serve as a defence and warning mechanism for the body.
- Do not apply stimulation while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk of injury.
- If you are in the care of a physician, consult with your physician before using this device.
- If you have had medical or physical treatment for your pain, consult with your physician before using this device.
- If your pain does not improve, becomes more than mild, or continues for more than five days, stop using the device and consult with your physician.
- Do not apply stimulation over your neck because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Do not apply stimulation across your chest because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

WARNINGS

- Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Do not apply stimulation when in the bath or shower.
- Do not apply stimulation while sleeping.
- Do not use the device on children, if it has not been evaluated for pediatric use.
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Apply stimulation only to normal, intact, clean, healthy skin.
- In the application you need to be careful when in the immediate vicinity a mobile phone is switched on, as this may affect the output of the device.
- If you are pregnant, you should not place any electrodes directly over the uterus or on the belly. The reason is that the heart of the foetus may be affected (so far there are no medical reports of damage).
- Simultaneous connection of a PATIENT to a high frequency surgical device may result in burns at the site of the STIMULATOR electrodes and possible damage to the STIMULATOR.
- Operation in close proximity to a shortwave or microwave therapy device may produce instability in the STIMULATOR output.
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.

Warnings when using a Vaginal / anal probe

- Pregnant women should not be treated in the area of the uterus.
- People with extra-urethral incontinence (fistula, ectopic urethra) may not be treated with a probe.
- People with over-sized incontinence should not be treated with a probe due to a closure.
- People with severe urinary retention in the upper urinary tract should not be treated with a probe.

CAUTIONS / CONTRAINDICATIONS

- People with complete peripheral denervation of the pelvic floor should not be treated with a probe due to a closure.
- People with total or partial prolapsed uterus / vagina should be treated with extreme caution.
- People with urinary tract infections need to contact their doctor before starting with the treatment.
- If skin irritation develops, treatment should be temporarily interrupted. If problems persist, consult your doctor. Hypersensitivity may occur in individual cases. The problem usually disappears when probes or gel is exchanged for a different type.

CAUTIONS / CONTRAINDICATIONS

Precautions

- TENS is not a substitute for pain medications and other pain management therapies.
- TENS devices have no curative value.
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- The success of the therapy depends to a considerable extent on the patient's choice of physician and whether they are qualified for the treatment of patients suffering from pain.
- The long-term effects of electrical stimulation are unknown.
- You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium.
- Use caution if you have a tendency to bleed internally, such as following an injury or fracture.
- Consult with your physician prior to using the device after a recent surgical procedure, because stimulation may disrupt the healing process.


- Use caution if stimulation is applied over the menstruating or pregnant uterus.
- Use caution if stimulation is applied over areas of skin that lack normal sensation.
- Use this device only with the leads, electrodes, and accessories recommended by the manufacturer.

Contraindications

- Long-term use can occasionally lead to skin irritation in the area where the electrodes were placed.
- The effectiveness of the treatment depends greatly upon the patient being treated by someone who is qualified in dealing with patients suffering from pain.
- Potential negative reactions include skin irritation and burns caused by the electrodes.
- You may experience skin irritation and burns beneath the stimulation electrodes applied to your skin.
- You may experience headache and other painful sensations during or following the application of electrical stimulation near your eyes and to your head and face.
- You should stop using the device and should consult with your physician if you experience adverse reactions from the device.

Before you start:


- Make sure the batteries are installed correctly.
- Connect the cable to the electrode pads and then connect to the unit.
- Attach the electrode pads to the area of pain.
- Do not use the electrode pads if it is scratched or damaged in any way.

Check/Replace the battery

1. Remove the Belt Clip (**K**) by sliding it down the unit.
2. Open the battery compartment cover (**L**).
3. Insert 4 batteries (type: AAA) into the battery compartment. Make sure you are installing the batteries correctly. Be sure to match the positive and negative ends of the batteries to the markings in the battery compartment of the device.
4. Recover the battery compartment (**L**) cover if you want to use it.
5. Bring on the Belt Clip (**K**) back by sliding it upwards the unit.

**Caution:**

1. Remove the batteries if the device is not in use for long periods of time.
2. Do not mix old and new batteries or different types of batteries.
3. Warning: If batteries leak and come into contact with the skin or eyes, wash immediately with copious amounts of water.
4. Batteries must be handled by an adult. Keep batteries out of the reach of children.
5. Only batteries of the same or equivalent type are recommended.
6. Do not use rechargeable batteries.
7. Remove exhausted batteries from the unit.
8. Dispose of batteries safely according to battery manufacturer's instructions.

Connect the electrode pads to the cables:

Hold the cable plug and insert into electrode pad connections. Make sure there are no bare metal of the pins exposed.

**Caution:**

Always use the cable that supplied with manufacturer or distributor and use the electrode pads with CE mark, or are legally marketed in the US under 510 (k) procedure.

Connect cable to unit

Before proceeding to this step, ensure the device is completely switched off. Hold the cable connector and plug it into the jack (A/B) at the top of the unit.

**Caution:**


Do not insert the plug of the patient lead wire into any AC power supply socket.

Place electrode pads on skin

Attach the self-adhesive electrodes on your body recommended by your doctor / therapist. Before applying electrodes, be sure the skin surface over which electrode pads are placed is thoroughly clean and dry. Make sure the electrode pads are pressed firmly to the skin and make good contact between the skin and the electrode pads. Place the electrode pads over the skin, attach them properly, firmly, and evenly.

**Caution:**

1. Before applying the electrode pads to the body, ensure the skin surface is clean and free from lotions or moisturizers.
2. Do not switch the device on when self-adhesive electrode pads are not positioned on the body.
3. Never remove the pads from the skin whilst the device is still switched on.
4. For hygiene reasons we recommend that the electrode pads are replaced every 30 days.
5. Please make sure the device is completely turned off before connecting.
6. It is recommended that, at minimum, 40mm x 40mm self-adhering based, square electrodes are used at the treatment area.
7. For hygienic reasons, each patient should use their own set of electrodes.

Attach the vaginal / anal probe to the cables:

Hold the cable connector and feed it into the sockets on the probe. Make sure that the bare metal of the pins is not exposed.

Caution: Always use probes from **Promed**.


Connect the cable to the device

Before proceeding with this step, make sure that the device is turned off completely. Hold the cable connector and plug it into the jack (A/B) at the top of the unit.

**Caution:**

Insert the plug of the patient lead cable not into an AC outlet.

Insert the vaginal / anal probe:

Take the probe out of the box and clean it under running water. Apply a thin layer of lubricant on the probe, and run the probe slowly and gently into the vagina or the anus, indicating to the flange at the bottom of the probe between the labia, or at the orbicular muscle. Make sure while inserting the vaginal probe that the silver contact surfaces are on the left and right.

Caution:

1. Do not turn the device on when the probe is not positioned in the body.
2. Never remove the probe when the device is switched on.
3. The probe must be properly cleaned before and after each use due to the hygiene reasons.
4. Make sure that the device is completely switched off before connecting.
5. For hygiene reasons, each patient should use a separate probe.

Switch on the device:

To switch on the device, press On/Off switch **"ON/OFF" (C)** on the upper right side of the device. The following screen appears:

There are 5 pre-set and a free programmable user program. Press the **"P"**-button **(G)**, to select the programme that you want to use and the screen on the LCD-Display **(F)** will look like this:

Selecting the programme mode

Each programme has 2 – 5 programme mode that you may select with **"S"**-button **(H)**.

Adjusting the channel intensity

The **Promed IT-6** has 2 output channels. You may either use a channel having a pair of electrodes, or both channels with two pairs of electrodes. There is an intensity regulator at the side of the device for each channel (**D**) and (**E**). By rotating each slider upwards, the intensity is increased when the intensity reaches to 1 mA, the device starts to operate on that channel. The maximum output voltage of the device is 90 mA. By rotating each slider downwards, the intensity is decreased when the intensity reaches to 0 mA, the device stops to operate on that channel.

Caution:

1. The intensity of the stimulation depends on the individual need of the patient.
2. The words „Work“ or „Rest“ on LCD-Display means that the device gives a signal or there is no signal.
3. If the electrodes are not placed firmly on the skin, or the device is not connected to the electrodes, the intensity automatically resets to 0 mA, provided the intensity is at least 10 mA.
4. Before removing the electrodes, make sure that the device is switched off.
5. If the treatment time is up, the intensity is automatically adjusted to 0 mA.
6. In case of an emergency please press the **“P”**-button (**G**) to cancel the stop the treatment. The device will remain in pause mode. Press the **“S”**-button (**H**) to continue the treatment.
7. It is recommended not to exceed a current density of 2 mA /cm², otherwise skin irritations or burns can occur. We recommend that:
 - Please use specific size electrodes provided by manufacturer to avoid

skin irritations or burns.

- For smaller electrodes, the maximum current setting of the waveform should be appropriately reduced.

- Please place the electrodes carefully, ensure that the entire surface of the electrode has good contact with the skin.

Programming the user program

Press **“ON / OFF”** (**C**) to switch on the device. Then press the **“P”**-button (**G**), to select the user programme PC1, it consists of 3 programme modes. You can set the parameters (pulse rate, pause time, pulse width, treatment time, and working time, ramp up time, ramp return time, sleep time, and waveform) according to your needs for each program mode.

Press the **“S”**-Button (**H**), to toggle between the individual program modes. The display of the programme phase on LCD-Display (**F**) will appear as follows:


Press and hold the **“S”**-button (**H**) to reach to the setting mode. Then either press the **“▲”**-button (**I**) or the **“▼”**-button (**J**), to select „Work Rest“ or „Cont“. If you select „Cont“, the stimulation is continuous and you can set the parameters like pulse frequency, pulse width and treatment time. If you select „Work Rest“, you can set all the parameters.


Press the **"S"**-button (**H**) and the wave form display flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to select the desired wave form (SY or AL).


Press again the **"S"**-button (**H**) and the pulse-width display flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to the desired pulse width from 50 μ s to 450 μ s.

Press again the **"S"**-button (**H**) and the pulse-display flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to set the desired pulse frequency from 3 Hz to 100 Hz.


Press again the **"S"**-button (**H**) and the treatment time display of the programme mode flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to set the desired treatment time of the programme mode from 1 min to 60 min.

Press again the **"S"**-button (**H**) and the pulse frequency-display flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to set the desired relaxation pulse frequency from 0 Hz to 4 Hz.


Press again the **"S"**-button (**H**) and the working hours-display flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to set the desired working hours from 2 seconds to 99 seconds.


- DE
- EN
- FR
- IT
- ES
- NL
- RU
- PL
- SE
- FI


Press again the **"S"**-button (**H**) and the ramp-up time display flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to set the desired ramp-up time from 1.0 seconds to 9.0 seconds.

Press again the **"S"**-button (**H**) and the ramp-return time display flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to set the desired ramp-return time from 1.0 seconds to 9.0 seconds.

Press again the **"S"**-button (**H**) and the rest time display flashes. Now either press the **"▲"**-button (**I**) or the **"▼"**-button (**J**) to set the desired rest time from 2 seconds to 99 seconds.

When you are done with the settings of the program mode, press the **"P"**-button (**G**), to confirm or save the settings. The device now switches to the waiting-mode. Now press the **"S"**-button (**H**) to another programme mode. Repeat all previous steps to set the parameters. When you're done with all the settings you can start the programme by regulating the intensity.

Functions of the back light

There are 3 colours of the back light that indicates the different states:

- Back light green: Normal state
- Back light blue: The intensity of the used channel is greater than 65 mA. With the further increase in the intensity, you must to exercise extreme caution.
- Back light red: The electrodes are not fixed to the skin or the device is not connected to the electrodes, and the set intensity is greater than 10 mA. The intensity is thus automatically reset to 0 mA, and the back light keeps flashing in red.

Checking the memory status

Press and hold the **"P"**-button (**G**) to reach to the saving mode. You can now read the entire treatment time and average intensity of channel 1 or 2 with the **"▲"**-button (**I**) or the **"▼"**-button (**J**). Press the **"P"**-button (**G**) again to return to waiting mode.


Caution:

If the total treatment time is less than 60 minutes, the unit will show you this.

At a total treatment time of up to 60 minutes, the total treatment time and average intensity is displayed.

The total treatment time is displayed up to a maximum value of 99:59 minutes.

To delete the memory, hold the **"P"**-button (**G**), pressed in display, the display „CLR“ flashes.

If you wish to delete all saved data, please press the **"S"**-button (**H**), to confirm the action. If you do not wish to delete the data, press the **"P"**-button (**G**) and the device switches back to waiting mode.

Safety lock

The security lock is activated automatically if no button is pressed for more than 30 seconds. To activate the safety lock itself, press the **"S"**-button (**H**) and the **"▲"**-button (**I**) simultaneously. This is a safety feature to prevent accidental changes of the settings, especially to prevent the set intensity. To deactivate the safety lock, press the **"S"**-button (**H**) and the **"▲"**-button (**I**) again. By pressing the **"S"**-button (**H**) during the treatment, you reach the pause mode and by pressing **"P"**-button (**G**) button during the treatment, you can stop the treatment.

Battery status display

If the battery status indicator  appears on the display, the used batteries must be replaced with new ones as soon as possible. The device can of course be operated for more hours with the illuminating display.

Disconnecting the device

Push the switch **"ON / OFF"** (**C**) down, to switch off the device.

Caution:

1. If you see the device in waiting mode and no button is pressed for more than 3 minutes, the display disconnects automatically to save the power.
To indicate the user thereon, that he has not pressed the button for long, a beep sound is heard after one minute, after second minute, two beeps and after 3 minutes, 3 beeps and the display is turned off.
By pressing any button the display is switched on again.
2. You must turn off the device to restore it to the factory settings. Now first press the button (**J**) and then push the **"ON / OFF"** (**C**) switch upwards.

Therapy-programme

There are 5 preset programmes and a fully programmable user program. In each program, up to 5 program modes can be adjusted. Please see the details in the table below:

Preset programmes

Programme	Programme mode	Wave form	Treatment time mode (min)	Frequency (Hz)	Pulse width (μ s)	Total treatment time (min)	Applies in the following incontinence type
P1	1	Synchronous (SY)	8	35	220	28	Urge incontinence
	2		20	10	240		
P2	1	Alternating(AL)	5	3	250	27	Mixed /faecal incontinence
	2		6	10	220		
	3		6	20	220		
	4		5	35	200		
	5		5	10	220		
P3	1	Pulse width modulation	5	20	300-450	15	Stress incontinence
	2	Modulation	5	3-35	175-275		
	3	Continuous	5	50	220		
P4	1	Continuous	8	10	240	25	Mixed /faecal incontinence
	2	Synchronous (SY)	12	20	275		
	3		5	35	200		
P5	1	Synchronous (SY)	20	10	240	38	Urge incontinence
	2		10	20	275		
	3		8	35	220		

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Free programmable user programme

Pro-gramme	Programme parameter	Adjustable values	Applies in the following incontinence type
PC1	Treatment time (min)	1 min~60min	All forms of incontinence depending on the set parameter
	Work mode	Work Rest/Cont.	
	Frequency	3 Hz~100 Hz	
	Rest frequency	0 Hz~4 Hz	
	Pulse width	50 μ s~450 μ s	
	Working hours	2 s~99 s	
	Ramp-up time	1 s~9 s	
	Ramp-down time	1 s~9 s	
	Rest time	2 s~99 s	
	Wave form	Synchronous (SY) / Alternating (AL)	

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Pay attention to the following before first time operation:

- The Promed IT-6 device is only intended for domestic use.
- Scrupulous supervision is required when the device is being used in the vicinity of children.
- Never place or use the unit in a wet or moist environment.
-  Do not use the unit under water, e.g., under the shower.
- Never place the unit in direct contact with fire, gas or oxygen, or hot objects such as hot plates.
- Take every precaution possible to ensure that the unit does not fall down or get damaged in any other way.
- If problems do arise with the unit, take it immediately to get it repaired.
- Do not lubricate or wash the unit.

 **Danger!**

- Do not place the unit in direct contact with water or other fluids.
- It is not permitted to use the unit outdoors.
- Never touch the unit with wet hands.
- Do not store the unit near a wash basin or bath tub as there is the danger that the unit may fall or be pulled into the wash basin or the bath tub.

 **Warning!**

- Never leave the unit unattended when children or inexperienced persons are in the vicinity of the unit.
- Make sure that children do not play with the device.
- Only use the unit for the applications described in this manual.
- This unit is not designed for use by people (including children) who have limited physical, sensory, or intellectual abilities or who do not have sufficient experience and/or sufficient knowledge. This does not apply if these people are supervised by a person responsible for their safety or receive instructions from them on how to use the unit.
- Do not use the unit underneath blankets or cushions as excessive heat may cause fire, injury or an electrical shock.
- Do not carry the unit by its cable and do not use the cable as a handle.
- It is not permitted to use the unit outdoors.
- When you have finished treatment, place the intensity controllers (**D und E**) in the "OFF" position.
- Never let children play with the packaging material; there is a risk of suffocation.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Problem	Possible Cause	Solution
Displays fail to light up	Battery contact failure	1. Try fresh batteries.
		2. Ensure batteries are inserted correctly. Check the following contacts: <ul style="list-style-type: none"> • All contacts are in place. • All contacts are not broken.
Stimulation weak or cannot feel any stimulation	Electrodes 1. Dried out or contaminated 2. Placement	1. Replace 2. Reattach the electrodes
	Lead wires Old/worn /damaged	Replace
	The intensity is too weak.	Use a higher intensity level.
Stimulation is uncomfortable.	Intensity is too high	Decrease intensity.
	Electrodes are too close together	Reposition the electrodes.
	Damaged or worn electrodes or lead wires	Replace
	Electrode active area size is too small.	Replace electrodes with ones that have an active area no less than 16.0 cm ² (4cm*4cm).
	Mayn't operate the device according to the manual.	Please check the manual before use.
Intermittent output	Lead wires	1. Verify connection is secure. Insure firmly.
		2. Turn down the intensity. Rotate lead wires in socket 90°. If still intermittent, replace lead wire.
		3. If still intermittent after replacing lead wire, a component may have failed. Call the repair department.

Stimulation is ineffective.	Misplacement of the electrode and probe	Repositioning electrode or probe. Contact clinician.
The skin becomes red and/or you feel a stabbing pain	Use the electrodes on the same site every time.	Re-position the electrodes. If at any time you feel pain or discomfort stop use immediately.
	The electrodes aren't stuck onto the skin properly.	Ensure the electrode is stuck securely on the skin.
	The electrodes are dirty.	Clean the electrode pads with a damp, lint free cloth or replace new electrode pads.
Output current stops during	The electrode/probe surface is scratched.	Replace probe/against a new electrode.
	The electrode pads come off the skin.	Turn off the device and stick the electrode pad firmly to the skin.
	The cable is disconnected.	Turn off the device and connect the cable
	The power of the batteries has been exhausted.	Please replace them with new batteries.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

- Do not place the unit in direct contact with water or other fluids.
- Switch the unit off when it is not in use.
- Before putting the **Promed IT-6** into storage for a longer period of time, remove the batteries from the compartment. Leaking batteries can damage the unit.
- Keep the unit and its accessories in a cool, dry place in the case provided.
- Handle the electrodes by their edges when removing them. To avoid damaging them, never pull at the electrode leads themselves.
- Do not make any sharp kinks in the connecting leads or electrodes.
- After use, stick the electrode pad onto the protective plastic film.
- Do not expose the device to direct sunlight and protect it against dirt and moisture.
- Never place any heavy objects on the device.
- You can clean your **Promed IT-6** by wiping it over carefully with a cloth moistened with mild soapy water. You may also use isopropyl alcohol or a soap solution. Household detergents and cleaning products are not suitable.
- Only use electrodes / probes approved by Promed for the device.
- A vaginal / anal probe is meant to be used only by one person. If used by more persons, then there is a high risk of cross contamination.
- Wash vaginal / anal probe in mild soapy water, rinse thoroughly with clean water and air dry with a clean lint-free cloth / tissue.
- Never put a not completely dried probe for storage in a plastic bag.
- The unit must not be returned to your dealer or the manufacturer for inspection and recalibration. If you need such inspections or recalibrations because of your internal Quality Management System don't hesitate to contact Promed.
- If you should experience other problems, consult your dealer, returning the unit if necessary. Never attempt to repair a fault yourself.

Power tools, accessories, and packaging should be recycled in an environmentally suitable manner.

Only for EU countries: Do not throw power tools in with the domestic waste! In accordance with the European Directive 2012/19/EU on waste electrical and electronic equipment and its implementation into national law, nonfunctioning  power tools must be collected separately and recycled in an environmentally suitable manner.

This symbol indicates that disposal of this product in the domestic waste is not permitted within the EU. Waste devices contain valuable recyclable material that should be recycled. Recycling also prevents uncontrolled waste disposal from damaging the environment and human health. Therefore, please dispose of all waste devices using the appropriate collection systems or send the device to the place of purchase for disposal. The place of purchase will then recycle the device.

Rechargeable batteries/batteries:

Do not throw your rechargeable batteries/batteries in the domestic waste, in fire, or in water. Rechargeable batteries/batteries should be collected, recycled, or disposed of in an environmentally friendly manner.

Only for EU countries:

According to the directive 91/157/EEC, defect or used rechargeable batteries/batteries must be recycled. Waste rechargeable batteries/batteries can be directly handed to:

Promed GmbH, Lindenweg 11, D-82490 Farchant.

This product was produced with the utmost care and was thoroughly tested before it left the factory. Therefore, we provide a warranty period of 24 months from the date of purchase for this product. **Promed** products conform to their description and the relevant specifications; it is your responsibility to ensure that the product which you have purchased is designed for your intended use.

- In the event of material defects and/or manufacturing errors which appear during proper use and are detected during the warranty period, we will replace all defective parts of the product free of charge within the warranty period, including the labour costs for warranty repairs.
- **The warranty does not cover:**
 - Normal wear on the product
 - Defects which were caused during transportation or storage of the product
 - Defects or damage which were caused through improper use or inadequate maintenance
 - Damage resulting from disregarding the notes in the operating manual
 - Damage resulting from modifications to the product which were not carried out by **Promed**.
 - Damage caused by sharp objects or as a result of torsion, compression, being dropped, an abnormal impact or other actions which are beyond **Promed**'s reasonable control.
 - Parts which are subject to wear (e.g. moving parts such as ball bearings, Motors, motor carbon brushes etc., fasteners, etc.) are generally excluded from the warranty.
- **A claim under the warranty is not enforceable if:**
 - the product is not sent back in its original packaging or in adequately secure packaging;
 - it has been modified or repaired by a person or company other than **Promed** or one of **Promed**'s authorised dealers;
 - the product has been repaired using parts which are not approved by **Promed**;
 - the serial number / LOT number of the product has been removed, deleted, changed or made illegible.
- For reasons of hygiene, any products which have been exposed to direct contact with a body or with bodily fluids (e.g. blood) must be packaged in an additional plastic bag before they are sent back. In this case, a specific note of this situation must be made in the package or in the accompanying documents.
- For components or products which are replaced during the warranty period, a warranty is only provided for the remainder of the original warranty period; provided that the product or part replacement was carried out by **Promed** or one of **Promed**'s authorised dealers.
- The warranty period begins on the day of purchase. Claims under the warranty must be established within the warranty period. Complaints arising after the expiry of the warranty period cannot be considered.
 - Under these warranty conditions, the warranty shall only apply if the date of purchase can be proven by means of a sales receipt or similar.
 - We reserve the right to make technical and visual changes as well as changes to the features.
- This warranty is only valid and enforceable in the country in which it was purchased by the original buyer, provided that it was **Promed**'s intention that the product should be offered for sale in this country. This warranty is also enforceable in any country within the European Economic Area in which **Promed** has an authorised importer or sales partner. Depending on each country, specific and differing warranties and guarantees may be applicable owing to the particular applicable legislation. These legal provisions are neither excluded nor limited by these warranty conditions. To the extent permissible under national law, the warranty period shall not be extended, renewed or otherwise affected by any subsequent resale, repair or replacement of the product.
- The provisions of the UN Convention on Contracts for the International Sale of Goods shall not apply.
- Our warranty conditions do not affect the seller's legal warranty obligation.
- To the fullest extent possible under the applicable mandatory legislation, these warranty conditions constitute your sole and exclusive legal remedy and shall apply in place of any other express or implied warranty conditions. **Promed** is not liable for any unusual, incidental, punitive or consequential damages including, but not limited to loss of profits, loss of use, loss of revenue, costs of replacement equipment or facilities, third party insurance claims, and damage to property which are the result of the purchase or use of the product or occur owing to a breach of warranty, a breach of contract, negligence, product defects or any other legal or statutory circumstances, even though **Promed** was aware of the possibility of such damages. **Promed** is not liable for any delay in warranty claims.
- **Promed** cannot be held liable for any possible translation errors.
- The following information is vital for problem-free processing:
 1. Original sales slip/receipt or dealer's stamp with the date of purchase
 2. Defect detected
 3. Unit name / type / serial / lot number

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EN | The requirements of the EC Directive 93/42/EEC in conjunction with EC Directive 2007/47/EC for class IIa medical products has been fulfilled according to the Standard EN 60601-1-2:2015.

CE according to the EC Directive 2014/30/EU and EC Directives 2011/65/EU, 1907/2006/EC, 2005/69/EC and 93/68/EEC.
The CE marking for the appliance refers to the EC Directive 93/42/EEC in conjunction with EC Directive 2007/47/EC.
Appliance classification: Class IIa.

Furthermore, the device fulfils the requirements of the standards:

EN 60601-1:2006 +A1:2013
EN 60601-1-2:2015
EN 60601-1-6:2010 +A1:2015
EN 60601-1-11:2015
EN 1041:2008 +A1:2013
EN ISO 15223-1:2016
EN ISO 10993-1:2009 +AC:2010
EN ISO 10993-5:2009
EN ISO 10993-10:2013
EN 62304:2015
EN 62366-1:2015

Electromagnetic disturbances are neither emitted from the appliance nor affect the appliance.

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Only transport in the original packaging.

Transportation/storage conditions

Temperature:
-10 °C – +55 °C
(+14 °F to +131 °F)

Relative humidity:
10 % - 90 % non-condensing
Atmospheric pressure:
700 to 1 060 hPa

Operating conditions:

Temperature:
+5 °C – +40 °C
(+41 °F to +104 °F)

Relative humidity:
30 % to 75 % non-condensing
Atmospheric pressure:
700 to 1 060 hPa

Service life of the device: 3 years
Service life of battery: new batteries will last for approx. Two week (when used for once a day, 45 mA intensity)
Service life of electrode pads: Pads can be cleaned and reused for up to 10~15 times.

An abrupt change in temperature is to be avoided as condensation may arise in the unit. Only operate the unit after the temperature has been equalized.

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- DE
- EN
- FR
- IT
- ES
- NL
- RU
- PL
- SE
- FI

- This device should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this device should be observed to verify normal operation in the configuration in which it will be used.
- Use of accessories other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- When the operating environment is relatively dry, strong electromagnetic interference usually occurs. At this time, the device may be affected as follows:
 - the device stops output;
 - the device turns off;
 - the device restarts;

The above phenomenon does not affect the basic safety and essential performance of the device, and the user can use it according to the instruction. If you want to avoid the above phenomenon, please use it according to the environment specified in the manual.

Table 1

Declaration - electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for domestic establishment and in establishment directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2

Declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Table 3

Declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3V 0.15 MHz to 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz	Not applicable	<p>Portable and mobile RF communications equipment should be used no closer to any part of device, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d=1.2\sqrt{P}$ 150 KHz to 80 MHz $d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

- Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which device is used exceeds the applicable RF compliance level above, device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating device.
- Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and device, as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	0.15 MHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.7 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Shenzhen Dongdixin Technology Co., Ltd.
Shanghai International Holding Corp. GmbH (Europe)
MedNet SWISS GmbH, Bäderstrasse 18, 5400 Baden,