
DOP1.1 DE | EN | ES | FR | SV | DA

GEBRAUCHSANWEISUNG | OPERATING INSTRUCTIONS | INSTRUCCIONES DE MANEJO
MODE D'EMPLOI | BRUKSANVISNING | BRUGSANVISNING

DEEP OSCILLATION® Personal

PHYSIOMED®
technology for therapy

DE

EN

ES

FR

SV

DA

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DA

Legal Notice

Distribution, reproduction and translation of the software and its documentation (or excerpts thereof) are prohibited without the prior written consent of PHYSIOMED ELEKTROMEDIZIN AG.

PHYSIOMED ELEKTROMEDIZIN AG reserves the right to change the software and associated data as well as the documentation without notice. All other rights reserved.

Manufacturer's Address

PHYSIOMED ELEKTROMEDIZIN AG
Hutweide 10
91220 Schnaittach
Germany
Phone: +49 9126 2587-0
Fax: +49 9126 2587-25
E-Mail: info@physiomed.de
Internet: http://www.physiomed.de

DEEP OSCILLATION® Personal is a registered trademark of PHYSIOMED ELEKTROMEDIZIN AG. The therapy mode has been patented worldwide.

DEEP OSCILLATION® Personal is made in Germany in compliance with the quality requirements of EN ISO 13485 and the applicable safety standards of the European regulation 2017/745 on medical devices.

The compliance with the regulations mentioned here is indicated by the CE marking on the device. The declaration of conformity might be requested from the manufacturer at the address given above.

The device has the CE marking 0123.

Operating Instructions last updated: 2023-11-10

Revision: 3.0

Table of Contents

Chapter 1	Instructions for Use	1
Chapter 2	Quick Reference Guide	2
2.1	Preparations, Start-up and Treatment.	2
Chapter 3	Introduction	3
3.1	Conventions Used	3
3.2	Device Overview.	4
3.3	Rating Plate	5
3.4	General Notes.	6
3.5	Essential Performance	6
3.6	Device Description.	7
3.7	Intended Purpose	7
3.8	Target Patient Group	7
3.9	Medical Indications and Claims.	7
3.10	Contraindications	8
3.11	Risks and Adverse Events.	8
3.12	Intended User	9
3.13	Intended Point of Use	9
3.14	Intended Duration and Intended Frequency of Use	9
3.15	Intended Location of Use and Intended Type and Duration of Contact	9
3.16	Precautions for Use	9
Chapter 4	Controls and Indicators	10
4.1	Display <1>	10
4.2	Data Selector <2>	11
4.3	Ports <3>	11
4.4	Card Reader <4>	11
4.5	Power Adapter Socket <5>.	11
4.6	Control Light <6>	12
Chapter 5	Device Operation	13
5.1	Battery Operation	13
5.1.1	Economy Mode	14
5.2	Preparations and Start-Up	14
5.3	Function Check	14

5.4	Hand Applicator	15
5.5	Pin Applicator	16
5.6	Therapy Card	16
Chapter 6	Treatment Preparation.	17
6.1	Device Errors	18
Chapter 7	Indications	19
Chapter 8	Individual Programs (only DEEP OSCILLATION Personal Pro)	21
Appendix A	Reference Information.	23
A.1	Customer Service & Maintenance.	23
A.2	Cleaning and Disinfection	23
A.3	Service Life and Disposal	23
A.4	Electromagnetic Compatibility	24
A.4.1	Environment Emission Limits	25
A.4.2	Immunity	25
A.4.3	AC connection for the supply input.	25
A.4.4	Test specifications for immunity and covers to HF wireless communication equipment	26
A.5	Training.	27
A.6	Settings.	27
A.7	Technical Data.	28
A.7.1	Environment Conditions	29
Appendix B	Scope of Delivery and Accessories.	30
B.1	Scope of Delivery	30
B.1.1	DEEP OSCILLATION Personal Basic	30
B.1.2	DEEP OSCILLATION Personal Pro	30
B.2	Available Accessories	31
B.3	CE Designation	32
Appendix C	Manufacturer's Recommendations	33
Appendix D	Warranty.	34
	Index.	35

Chapter 1 Instructions for Use

Important

Please carefully read the instructions and warnings as well as the contraindications before first using the device.

To prevent the risk of strangulation or suffocation the medical device and all supplied accessories must be kept out of reach of children.

Please contact PHYSIOMED ELEKTROMEDIZIN AG (refer to [Manufacturer's Address](#) on page 2) or consult a qualified person active in health services if you have problems when first using the device, or observe any unusual function during operation.

The service life of the device is limited to 8 years. Please observe the following instructions:

- [Scope of Delivery and Accessories](#) on page 30
- [Cleaning and Disinfection](#) on page 23

A fully charged battery typically provides capacity for up 24 treatments or 8 hours of operation. The typical service life of a battery is about 500 charging cycles.

Important

No treatment is possible while charging the battery.

How to reach the maximum service life of your battery

- (a) Charge the battery before it has been completely discharged.
- (b) Never charge a completely charged battery again. Overcharging reduces the service life of the battery.
- (c) Charge and store the battery at a room temperature of 20°C to 40°C.
- (d) Charge the battery when you haven't used the device for a longer timespan.
- (e) Fully charge the battery after every 25th charging cycle.

Chapter 2 Quick Reference Guide

With DEEP OSCILLATION® Personal you have acquired an extremely versatile deep oscillation system. The device will only show its true potential, however, if you are well informed about its functions. For this reason, carefully read these Operating Instructions and familiarize yourself with the use of the device.

You will find an overview of the symbols used in these instructions or applied on the device under [Conventions Used](#) on page 3 and [Rating Plate](#) on page 5.

The structure of the device and the function of the individual elements are described under [Device Overview](#) on page 4 and [Controls and Indicators](#) on page 10.

The device has the effects described under [Medical Indications and Claims](#) on page 7. Before treatment please observe [Contraindications](#) on page 8 .

You will find instructions on how to clean the device under [Cleaning and Disinfection](#) on page 23 .

The device has to be maintained regularly (refer to [Customer Service & Maintenance](#) on page 23).

2.1 Preparations, Start-up and Treatment

How to start up the device

- (1) Prepare treatment as described under [Treatment Preparation](#) on page 17.
- (2) Insert the therapy card with the chip side facing down into the **Card Reader <4>**.

NOTICE

If the card is defective or incompatible, a respective error message will be displayed.

- (3) Press the **Data Selector <2>** for about one second. You will hear a signal tone, the **Display <1>** will be illuminated and the welcome screen with the start menu will be displayed.
- (4) You can now select an indication and begin the treatment (refer to [How to perform a treatment with indications](#) on page 19).

How to switch off the device

- (1) Press the **Data Selector <2>** for about five seconds. The **Display <1>** goes out.
- (2) Remove the therapy card and unplug the cables from the **Ports <3>**.

Chapter 3 Introduction

3.1 Conventions Used

Please note the following typographic conventions used in these operating instructions:


- Cross references and important terms used for the first time in this document are written in *italics*.
- Names of menus and symbols on the display are written in **bold typeface**.

Paragraphs that deserve special attention are highlighted in the following way:

Symbol	Type	Meaning
	Tip	Gives extra hints for more convenient operation
	Note	Provides background information for better understanding
	Important	Prevents misunderstandings that might lead to limited operation of the device or insufficient therapeutical results
	Caution	Alerts in cases of possible damage to the device or risks of injury

3.2 Device Overview

Front


Legend

1	Display	4	Card Reader
2	Data Selector	5	Power Adapter Socket
3	Ports	6	Control Light

Bottom


Top

3.3 Rating Plate

The rating plate is on the rear side of the device:

The symbols on the rating plate have the following meaning:

Symbol	Meaning
	Manufacturer address with year and month of production
	This is a medical product
	Serial number: The first two digits of the number indicate the year and the last two the month of the date of production.
	Unique identifier of the medical product
	CE marking of the notified body 0123
	The device contains a Lithium ion battery. Do not expose it to excessive heat or open fire. The battery can only be replaced by PHYSIOMED ELEKTROMEDIZIN AG. Never open the device to replace the battery by yourself.
	Separate disposal of electronic and electric devices. Used appliances have to be disposed of separately from the domestic waste.
	Rating according to IEC/EN 60601-1: Application part type BF (body floating with applied current) An application part type BF is connected to the patient to transmit electrical energy or an electrophysiological signal to or coming from the body.
	Protected against solid particles (e.g. dust or dirt particles) with a diameter ≥ 12.5 mm; protection against falling drip water
	CAUTION! Follow the Operating Instructions and consider the physiological effects of the selected parameters.

Symbol	Meaning
	Follow operating instructions
	Do not charge during treatment! To charge, only use the supplied power adapter (Ref no. 55965).

3.4 General Notes

DEEP OSCILLATION® Personal complies with the technical standard IEC 60601-1, and is assigned to class IIa according to the Medical Devices Regulation (2017/745).

It is not intended for operation in explosion hazard zones or hydrotherapy rooms. Drastic temperature changes should be avoided, since condensation could be caused within the device. Do not start up the device until it is in thermal equilibrium with its environment!

The device may only be operated by intended users, and only in accordance with the Operating Instructions. Operating the device in the proximity (e.g. 1 m) of a short-wave or micro-wave therapy unit may cause output irregularities and should be avoided for this reason, as well as simultaneous connection of the patient to high-frequency surgical device.

The device must be stored in a way so that dust and particles cannot enter the ports. The device must be protected from falling down.

The device must not be stored or operated at elevated temperatures. Direct exposure to sunlight must be avoided. The maximum temperature for operation is 40°C. For every degree of excessive temperature, three minutes time for accommodation have to be assumed as a guideline.

CAUTION

Technical modifications to the device are not admissible and will inevitably lead to loss of approval as a medical device.

The assembly of medical electrical systems (ME systems) and their modification during their effective service life require verification of the compatibility with the standard for ME systems (EN60601-1, 3rd Edition) by the responsible organization (system builder).

Any serious incident that occurred when operating DEEP OSCILLATION® Personal has to be reported to PHYSIOMED ELEKTROMEDIZIN AG and to the responsible authorities of the respective country.

3.5 Essential Performance

The essential performance is a pulsed electrostatic field output for a massage/deep oscillation effect.

3.6 Device Description

DEEP OSCILLATION® Personal is a deep oscillation system which is successfully used in various areas of medicine.

DEEP OSCILLATION® Personal has two modes of operation:

- **Treatment:** In this mode, the device is separated from the mains. When the power adapter is plugged in, the device cannot be switched on and treatment can not be continued. Connecting the power adapter with the mains will reduce the intensity automatically to zero, treatment will be interrupted, and the device will be switched off.
- **Charging:** While charging, treatment is not possible.

CAUTION

During the charging procedure, the patient must neither be connected with the titanium neutral element, nor with the hand applicator!

DEEP OSCILLATION® Personal is available in two variants:

- DEEP OSCILLATION® Personal Basic
- DEEP OSCILLATION® Personal Pro allows to additionally select the treatment frequencies.

DEEP OSCILLATION® Personal had been registered by PHYSIOMED ELEKTROMEDIZIN AG to treat conditions in human medicine.

The functionality of DEEP OSCILLATION® Personal is based on a pulsed electrostatic field built up in the patient's body region. Its frequency changes between 5 - 250 Hz, according to the selected indication. Due to the movement of the hand applicator or special glove, a vibrating effect with deep impact is induced in the patient's tissue.

3.7 Intended Purpose

DEEP OSCILLATION® Personal is an instrument intended for delivering transcutaneous pulsed electrostatic therapy.

3.8 Target Patient Group

The target group of DEEP OSCILLATION® Personal are adult patients (18 years and above).

3.9 Medical Indications and Claims

Treatment with DEEP OSCILLATION® Personal is indicated for the following medical indications:

Indication	Benefit Claim
Acute musculoskeletal pain (ICD-10 WHO: FC0Z)	reduces acute musculoskeletal pain
Fibromyalgia (ICD-10 WHO: MG30.01)	reduces fibromyalgia pain

Indication	Benefit Claim
Lymphedema (ICD-10 WHO: BD93)	reduces lymphedema volume
Lipedema (ICD-10 WHO: EF02.2)	reduces lipedema volume
Chronical obstructive pulmonary disease (stable phase) (ICD-10 WHO: CA22)	improves lung function
Skin burns (max. second degree affecting less than 15% of body surface) (ICD-10 WHO: ND9Z)	promotes skin burn healing
Raynaud Phenomenon (ICD-10 WHO: BD42)	alleviates symptoms
Parkinson's disease (ICD-10 WHO: 8A00.0)	improves mobility
Hamstring muscle tightness (ICD-10 WHO: not reported)	improves mobility

**WARNING**

DEEP OSCILLATION® Personal is not indicated for the following indications for which the clinical evidence in the literature is limited or is insufficient to demonstrate clinical benefit: fibrocyst, bronchial asthma, postpartum swelling, muscle fatigue, osteoarthritis, gonarthrosis, arthrosis, cellulite, scars fibrotic adhesion, esthetics applications (ageing, dark circles, etc.)

3.10 Contraindications

Therapy with DEEP OSCILLATION® Personal is not to be used in the following cases:

- acute infections
- acute inflammations with pathogenic germs participation
- active tuberculosis
- acute venous diseases (untreated thromboses)
- untreated malignant processes
- erysipelas
- patients and therapists with cardiac pacemakers and other electronic implants, cardiac insufficiency, decompensated cardiac edema, and cardiac arrhythmia
- pregnancy
- hypersensitivity to electrostatic fields
- infectious skin diseases
- vertebrobasilar insufficiency (VBI)

3.11 Risks and Adverse Events

The following residual risks of DEEP OSCILLATION® Personal might occur:

Summary of the residual risks and adverse events for DEEP OSCILLATION Personal

Undesired Effect	Classification	Frequency
Erythema, red discoloration of the skin	IMDRF code: E171601, Med-DRA:10015150: Erythema	very common (> 10% patients)
Itching sensation, an irritating feeling that produces the desire to	IMDRF Code: E1708, Med-DRA:10023084: Itching	very common (> 10% patients)

Undesired Effect	Classification	Frequency
scratch		
Discomfort, a feeling of mental or physical uneasiness, pain or distress	IMDRF Code: E2311, Med-DRA:10013082: Discomfort	very common (> 10% patients)

3.12 Intended User

The intended users of DEEP OSCILLATION® Personal are health professionals and patients.

3.13 Intended Point of Use

The intended point of use is home or dry rooms designed for medical purposes (hospitals, physical therapy, and practice rooms with patient beds).

Avoid explosion hazard zones, moist areas (such as bathrooms) or hydrotherapy rooms.

3.14 Intended Duration and Intended Frequency of Use

The intended duration and frequency are a maximum of three sessions a day and a maximum cumulative daily use of three hours.

Alternate each therapy session with the gloves with a session using the applicator handhold.

3.15 Intended Location of Use and Intended Type and Duration of Contact

DEEP OSCILLATION® Personal is for external use only (direct contact with intact skin). The duration of contact is the same as the duration of use.

Place surgical drapes on wounds and burns prior to the treatment if sterile conditions must be maintained.

Avoid the following locations:

- mucous membranes
- the neck close to the thyroid gland
- the region directly above the heart

3.16 Precautions for Use

- All parts of the oscillator head must be dry before being used again.
- Switch off the device beforehand and disconnect it from the mains voltage when cleaning it.
- Clean and disinfect the oscillator head after each treatment.
- Membranes are for single patient use only.
- Check and replace membranes that are soiled, damaged, worn with use and no longer thick enough.
- Replace the membrane at least every 25 applications.


Chapter 4 Controls and Indicators

The design of DEEP OSCILLATION® Personal allows for easy operation, despite a variety of functions. Due to its small size, the device is very easily transportable. It has been designed for use outside of therapy rooms, and it runs by a rechargeable battery.

All controls and indicators are integrated into the casing, thus allowing for easy cleaning of the device's surface. The device's microprocessor monitors the safety-related components, prevents from erroneous operation and checks all functions of the device after switching it on.

The following section introduces the individual indicators and controls of DEEP OSCILLATION® Personal. The numbers in angle brackets refer to the [Device Overview](#) on page 4.

4.1 Display <1>

In the **Display <1>** the available options or therapy parameters are shown, depending on the context.

The battery symbol on the upper right corner shows the current charging condition of the batteries (also refer to [Battery Operation](#) on page 13).

A parameter can be selected using the **Data Selector <2>**.

Symbols in the Display

	Battery Indicator (refer to Battery Operation on page 13)
	Indications
	Individual Programs (only DEEP OSCILLATION® Personal Pro)
	Settings Menu (refer to Settings on page 27)
	Button Back

4.2 Data Selector <2>

The **Data Selector <2>** serves to select the therapy parameters shown in the **Display <1>** as well as to adjust the intensity.


You can switch between the individual options by rotating the selector and activate the respective function by pressing it.

Tip

Reduce the intensity to zero and press the data selector three times quickly to return to the welcome screen or first menu level.

4.3 Ports <3>

The **Ports <3>** serve to connect the supplied cables (Ref. no. 00261 and 00262, or optionally 00355 or 00356) for the hand applicator or pin applicator and the adhesive electrode or Titanium neutral element with the device.

Instructions on how to connect the cables can be found in section [How to prepare the applicator](#) on page 17.

4.4 Card Reader <4>

The **Card Reader <4>** serves to insert the therapy card which enables access to the user menus. This is also where the therapy parameters are stored. The therapy card has to be inserted *with the chip side facing down*.

4.5 Power Adapter Socket <5>

The **Power Adapter Socket <5>** can be found on the bottom of the device. The supplied power adapter (Ref. no. 55965, Sinpro HPU15-105) to charge the battery is plugged in here.

4.6 Control Light <6>

The **Control Light <6>** under the **Data Selector <2>** lights up in blue as soon as the intensity is increased. In case of low therapy frequency, the light flashes in sync with the pulses. When the electrostatic field is cut short, the light extinguishes even with non-zero intensity.

A functional check is conducted when switching on (refer to [Function Check](#) on page 14). The LED is illuminated in green light to show that the device is ready to use.

NOTICE

The **Control Light <6>** will light up in yellow light in case of an error. The according error code is shown in the **Display <1>**.

Chapter 5 Device Operation

5.1 Battery Operation

Important notes when handling the battery

The device contains a Lithium ion battery. Don't expose it to excessive heat, direct sunlight, or open fire.

Never toss batteries into the fire, there's danger of explosion!

Avoid contact with skin and eyes, if liquids should leak out due to battery damage. Rinse the affected areas immediately with water and look for medical attendance.

Immediately contact our service when batteries had been damaged.

The battery can only be replaced by an authorized service partner. Never open the device to replace the battery.

Only use the supplied power supply (Ref. no. 55965, refer to [Scope of Delivery](#) on page 30) for charging and observe the instructions to charge the battery.

DEEP OSCILLATION® Personal uses rechargeable Lithium ion batteries for power supply, which allow for operation independently from the mains.

Important

In order to ensure a long service life, the battery must be charged completely when first charged. The first charging procedure should not be interrupted.

The charging status of the batteries is displayed on the **Display <1>**:

Battery completely charged

How to Charge the Battery

CAUTION

During the charging procedure, the patient must neither be connected with the titanium neutral element, nor with the hand applicator!

- (1) Plug the supplied power adapter into the **Power Adapter Socket <5>** on the bottom of the device. The battery symbol will start to flash to indicate that the battery is being loaded.
- (2) Unplug the power adapter from the device when the charging procedure is finished.

5.1.1 Economy Mode

The device automatically switches over to economy mode to save power. There are two stages:

- (a) After about one minute the **Display <1>** will be dimmed. Press the **Data Selector <2>** to reset the illumination to the value configured in the settings menu (refer to [Settings](#) on page 27).
- (b) After about five minutes the device will switch off completely. Start it up again as described under [Preparations and Start-Up](#) on page 14.

5.2 Preparations and Start-Up

How to Start up the Device

- (1) Insert the therapy card with the chip side facing down into the **Card Reader <4>**.

NOTICE

If the card is defective or incompatible, a respective error message will be displayed.

- (2) Press the **Data Selector <2>** for about one second. You will hear a signal tone, the **Display <1>** will be illuminated and the welcome screen with the start menu will be displayed.
You can now select an indication and begin the treatment (refer to [How to perform a treatment with indications](#) on page 19).

How to Switch off the Device

- (1) Press the **Data Selector <2>** for about five seconds. The **Display <1>** goes out.
- (2) Remove the therapy card and unplug the cables from the **Ports <3>**.

5.3 Function Check


A function check is always necessary when you are unsure about whether the device is working properly.

How to Check the Function of the Device

- (1) Press the **Data Selector <2>** for several seconds until the device switches off.
- (2) Press the **Data Selector <2>** again until the device switches on:
 - If the device is working properly, the available options are displayed on the **Display <1>**. Additionally the **Control Light <6>** turns green (refer to [Control Light <6>](#) on page 12). You can operate the device.
 - If an error code is displayed, there is a device error. Proceed as described under [Device Errors](#) on page 18.

5.4 Hand Applicator

The hand applicator consists of the following parts:

(1) Applicator handhold (Ref. No. 00379)

(2) Cover

(3) Oscillator Plate

(4) Membrane

(5) Clamp Ring

NOTICE

(2) - (5) make up the oscillator head in the respective size.

Hand Applicator (Exploded View)

The oscillator head is screwed into the handle. To remove the head from the applicator handhold, you have to rotate it to the left.

Important

Any damage to the oscillator membrane might lead to harmless (but unpleasant) sensations on the skin.

Replace the membrane immediately if damaged.

Replace the membrane after each person to be massaged.

How to Replace the Membrane

- (1) Take off the clamping ring and remove the membrane.
- (2) Insert the new membrane with the foil facing down loosely into the clamping ring, so that the membrane foil protrudes laterally.

Changing the Membrane - Inserting

- (3) Press the oscillator plate firmly onto the foam rubber area of the membrane until the clamping ring engages with the oscillator plate. Make sure that the foil is tightly fixed with the clamping ring on the outside, and that the clamping ring is engaged!

Changing the Membrane - Engaging the Oscillator Plate

5.5 Pin Applicator

The pin applicator (Ref. No. 00396) can be used in the same way as the hand applicator.

Pin Applicator

5.6 Therapy Card

The therapy card must be inserted to perform a treatment or to access the basic settings.

The therapy cards feature the following functions:

- Basic: Indications
- Pro: Indications and individual
- DEEP OSCILLATION: Memory card with one treatment (created using DEEP OSCILLATION-Evident)

CAUTION

The treatment card must not be copied, and any attempt to reproduce it or read out the data will invalidate it! PHYSIOMED ELEKTROMEDIZIN AG accepts no responsibility for damage incurred in this way!

Chapter 6 Treatment Preparation

CAUTION

Before each treatment, adjust the intensity to match the respective conditions and pay attention to the contraindications (refer to [Contraindications](#) on page 8)!

During the charging procedure, the patient must neither be connected with the titanium neutral element, nor with the hand applicator!

Important

You have to choose the intensity according to the minimum dose principle, i.e. in the lower effective range!

The hand applicator or special glove should stay in contact with the patient's skin and should be moved on the body surface without lifting or stopping. Distinguish between therapy movement and leading movement: The therapy movement (to the center of the body) has to be performed with the appropriate pressure, the leading movement (back from the center of the body) without pressure.

When treating with low frequencies (10 - 40 Hz), perform the therapy movement especially slow. The leading movement may be quicker.

After switching to a new therapy phase, it might be necessary to re-adjust the intensity to reach the same strength of oscillation!

Important

If the skin in the zone of contact reddens severely, treatment should immediately be stopped and the applicator should be removed from the skin. Please wait until the condition of the skin returns to normal. A rapidly fading away skin reddening below the electrode is harmless and can be explained by the locally increased blood circulation. If the skin irritation persists, however, and is possibly accompanied by itching or inflammation, you have to consult your physician before any further treatment.

NOTICE

The device operates in biphasic mode, therefore it is not necessary to care about the polarity when plugging in the cables.

How to prepare the applicator

- (1) Plug the *connection cables DEEP OSCILLATION®* (Ref. no. 00261) for the hand applicator and for the Titanium neutral element into the **Ports <3>**.
- (2) Connect the Titanium neutral element with the other end of the cable and let the patient grip it. Al-

ternatively, you can attach it on any part of the patient's body.

- (3) Connect the applicator to the appropriate cable.
- (4) Continue with the steps described under [How to perform a treatment with indications](#) on page 19.

How to prepare the treatment with special gloves

- (1) Plug the supplied *Connection cable grey for adhesive electrodes* (Ref. no. 00262) as well as one of both *connection cables DEEP OSCILLATION®* (Ref. no. 00261) into the **Ports <3>**.
For spiral leads: Plug the optional *Spiral lead for adhesive electrodes* (Ref. No. 00355) as well as one *Spiral lead universal* (Ref. No. 00356) into the **Ports <3>**.
- (2) Connect one *PHYSIOPADS adhesive electrode* (Ref. no. 00589) with the fastener button of the grey connection cable. Attach the adhesive electrode to the arm or leg of the therapist. Avoid sensitive skin areas and remove grease and skin residues to ensure good conductivity.
- (3) Connect the connection cable DEEP OSCILLATION® with the Titanium neutral element and have the patient grip it.
- (4) Pull the supplied *special gloves* (Ref. no. 00346) over the hands of the therapist.
- (5) Continue with the steps described under [How to perform a treatment with indications](#) on page 19.

6.1 Device Errors

If a functional error is detected during the automatic selftest routine or operation, one of the following error codes is shown on the **Display <1>**:

- **ERROR_01 ... ERROR_05**

In this case, treatment is automatically interrupted. Switch off the device by pressing the **Data Selector <2>**. If the error code persists even after several selftests, consult the technical service.


In case the device does not work properly, or the treatment stops with a 10 seconds continuous tone, please charge the batteries (refer to [How to Charge the Battery](#) on page 13). Your DEEP OSCILLATION® Personal should work properly afterwards. But if the treatment still stops with a continuous tone of 10 seconds, the device does not work properly, you notice some unusual operation, or any incident occurs, please also contact our technical service.

Chapter 7 Indications

With DEEP OSCILLATION® Personal you can perform treatments according to preset indications which are stored together with the corresponding therapy parameters on the therapy card. For each displayed indication you can refer to the recommended therapy instructions. Each indication features an image of a body region to be treated, which shows the preferable direction of treatment.


How to perform a treatment with indications

- (1) Before treatment, prepare the patient's skin area to be treated by drying or powdering.
- (2) Start up the device as described under [How to Start up the Device](#) on page 14.
- (3) Select **INDICAT** in the side menu to go to the list of indications. The menu shows the currently selected page (e.g. **1/9**).
- (4) Rotate the **Data Selector <2>** until you get to the page with the desired indication.
- (5) Press the **Data Selector <2>** to access the page.
- (6) Rotate the **Data Selector <2>** to navigate to the desired indication.
- (7) Press the **Data Selector <2>** again. The therapy parameters of the selected indication are displayed. The menu shows the currently selected page (e.g. **1/4**)

Indications - treatment parameters

- (8) Rotate the **Data Selector <2>** again. An image showing a body area to be treated as well as the treatment direction is displayed.

Indications - direction of treatment

- (9) On the subsequent page, instructions to conduct the treatment are shown.
- (10) Rotate the **Data Selector <2>** to return to the first page.
- (11) Press the **Data Selector <2>** two times so that the intensity is shown underlined in blue.
- (12) Rotate the **Data Selector <2>** to select the desired intensity. The current treatment phase is indic-

ated by a black dot. The treatment begins, the treatment time elapses and the **Control Light <6>** is illuminated.

In the example shown here treatment with two different therapy phases is performed. Perform the treatment and pay attention to the therapy instructions shown before.

In the case of treatments with more than one phase, therapy parameters are automatically switched after the end of the first therapy phase. As soon as the time for the last treatment phase has elapsed, the intensity decreases to zero. You can now perform one more treatment or switch off the device.

If the patient senses discomfort, the therapy duration shall be cut short or even suspended. Proceed in the following way:

- (1) Click the therapy time with the **Data Selector <2>** so that it is underlined.
- (2) Rotate the **Data Selector<2>** to the left and cut short the therapy time to the desired interval.
- (3) Press the **Data Selector <2>** to confirm the selected therapy time.

**Tip**

Press the data selector three times quickly to return to the welcome screen or first menu level. The intensity has to be zero to do so.

Chapter 8 Individual Programs (only DEEP OSCILLATION Personal Pro)

Individual programs allow you to adjust therapy parameters such as duty cycle and frequency. The corresponding **INDIVID** menu is only available for DEEP OSCILLATION® Personal Pro, and cannot be selected in the DEEP OSCILLATION® Personal Basic version.

You have two options for individual treatment:

- INDIVID Sequence
- INDIVID Manual

How to perform an INDIVID SEQUENCE treatment

- (1) Choose **INDIVID** in the side menu to get to the individual programs:


Individual Program - Page 1 - Selection of Parameters

- (2) Click again on **INDIVID** to get to the first page.
- (3) Click again on **INDIVID**. You can now rotate the **Data Selector <2>** to select the treatment duration, the duty cycle (e.g. **3:1**), the frequency (e.g. **200 Hz**) and the treatment duration of the three phases. Move the cursor, which is indicated by a horizontal bar, to the desired parameter. To select the value of a parameter, first press the **Data Selector <2>** to activate the parameter. Then rotate the **Data Selector <2>** to select the desired value.
- (4) When you increase the intensity, the currently active phase is indicated by a black dot. The sequence is completed phase after phase. Once the intensity is up you can only change the duty cycle and the frequency. The duration within a sequence cannot be changed. Still, you can abbreviate (but not extend) the duration of the ongoing treatment.
- (5) Click on the symbol < in the top row of the **Display <1>** to return to the overview page.

How to perform an INDIVID MANUAL treatment

- (1) Choose **INDIVID** in the lateral menu to get to the programs:

Individual Program - Page 2 - selection of parameters

- (2) Click on **INDIVID** once more to get to the second page.
- (3) Click on **INDIVID** again. You can now rotate the **Data Selector <2>** to set the intensity, the treatment duration and the duty cycle (e.g. **3:1**), the lower frequency (e.g. **50 Hz**), the upper frequency (e.g. **125 Hz**) and the duration of the three phases. Move the cursor, which is indicated by a horizontal bar, to the desired parameter.
To set the value of a parameter, first press the **Data Selector <2>** to activate the parameter. Then rotate the **Data Selector <2>** to set the desired value.
- (4) When you increase the intensity, the currently active phase is indicated by a black dot. During treatment you can change all of the parameters (duty cycle, upper and lower frequency, intensity and duration).
- (5) Click the symbol **<** in the top row of the **Display <1>** to return to the overview page.

Appendix A Reference Information

A.1 Customer Service & Maintenance

**CAUTION**

The device must never undergo service or maintenance during operation.

The manufacturer is responsible only for the original security properties of the device. The device must be used in accordance with the operating instructions.

The device must exclusively be maintained by service personnel authorised by PHYSIOMED ELEKTROMEDIZIN AG. A certification must be obtained from the service personnel about the type and extent of the work, possibly with information relating to the modifications in the rated data or the work range. The certification must also contain the date of execution as well as company name and signature. Defective components which influence the safety of the device must be replaced by original parts. Parts lists and service instructions can be made available to qualified service personnel of the user upon request.

For any questions regarding technical description please refer to the service personnel at PHYSIOMED ELEKTROMEDIZIN AG.

The device and its accessories must regularly be serviced. The defined safety-relevant inspection dates are listed in the manufacturer's recommendations (see [Manufacturer's Recommendations](#) on page 33).

A.2 Cleaning and Disinfection

Regularly clean both device and accessories using a surface disinfectant, e.g. 70 percent medical isopropyl alcohol or Dismozon® pur. Ensure that the disinfectant does not contain benzol, ethanol, alcohol (except isopropyl alcohol), organic agents or thinners. Prior to cleaning, always switch off the device and separate it from the mains.

Use a soft sponge cloth for cleaning. Take care that no liquid substances get into the device.

Clean and disinfect the oscillator head after each treatment. Check your accessories regularly. Replace the oscillation head not later than after 25 treatments.

**Important**

All parts of the oscillator head must be absolutely dry before using them next time!

A.3 Service Life and Disposal

Due to legislation, the service life of this medical device is limited to a maximum of 8 years. The membranes have a service life of 2 years.

The device has to be disposed of in compliance with legal regulation regarding electronic devices. It is not allowed to put used devices into your domestic waste.

Environmentally responsible disposal of devices when service life has elapsed can be done via the manufacturer.

A.4 Electromagnetic Compatibility

Declaration

Medical electrical equipment is subject to special precautions with regard to EMC and must be installed and commissioned in accordance with the EMC instructions contained in the accompanying documents. In particular, medical electrical equipment can be affected by portable and mobile HF communications equipment. The manufacturer guarantees the conformity of the device with the EMC requirements only if the accessories listed in the EC Declaration of Conformity are used. The use of other accessories may result in increased electromagnetic interference or a reduced resistance against electromagnetic interference.

In order to ensure the basic safety and the essential performance characteristics of the device with regard to the electromagnetic disturbances during the expected period of use, maintenance must be performed within the specified interval and according to the instructions in the technical service manual.

The device must not be placed directly adjacent to or stacked with other devices. If such an arrangement is nevertheless required, the device must be monitored to verify its intended operation in this arrangement. Further EMC information can be found under [General Notes](#) on page 6 as well as in the technical information on the following pages.

If there is a loss or a reduction of the essential performance characteristics due to electromagnetic disturbances, it is to be expected that the device will issue an error message in the event of more severe interference. The source of the error must be eliminated and the device be re-started.

According to the EMC regulations for medical devices, the manufacturer is required by law to provide the following information.

CAUTION

The use of accessories, converters and cables other than those specified may result in an increased emission or reduced immunity to interference.

Portable HF communications equipment including its accessories such as antenna cables and external antennas should not be used closer than 30 cm (or 12 inches) to the parts and cables of the device designated by the manufacturer. Non-observance can lead to a reduction in the performance characteristics of the device.

DEEP OSCILLATION® Personal has been tested for usage in professional health care facilities and areas of domestic health care.

A.4.1 Environment Emission Limits

The device is intended for operation in an environment as specified below. The customer or the user of the device should ensure that it is operated in such an environment.

Phenomenon	Professional health care facility	Area of domestic health care
Power-induced and radiated interference emissions	Class B, Group 1 (according to CISPR 11)	Class B, Group 1 (according to CISPR 11)
Distortions caused by harmonics	Refer to IEC 61000-3-2	Refer to IEC 61000-3-2
Voltage fluctuations and flicker	Refer to IEC 61000-3-3	Refer to IEC 61000-3-3

A.4.2 Immunity

The device is intended for operation in an environment as specified below. The customer or the user of the device should ensure that it is operated in such an environment.

Phenomenon	EMC basic standard or test method	Immunity test level for professional health care facility	Immunity test level for areas of domestic health care
Discharge of static electricity and connection to patients	IEC 61000-4-2	± 8 kV contact discharge ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air discharge	± 8 kV contact discharge ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air discharge
High-frequency electromagnetic fields	IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz
High-frequency electromagnetic fields in the immediate vicinity of wireless communication devices	IEC 61000-4-3	Refer to Test specifications for immunity and covers to HF wireless communication equipment on page 26	Refer to Test specifications for immunity and covers to HF wireless communication equipment on page 26
Magnetic fields with energy-technical nominal frequencies	IEC 61000-4-3	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz

A.4.3 AC connection for the supply input

The device is intended for operation in an environment as specified below. The customer or the user of the device should ensure that it is operated in such an environment.

Phenomenon	EMC basic standard	Immunity test level for professional health care facility	Immunity test level for areas of domestic health care
Fast transient electrical disturbances/bursts	IEC 61000-4-4	± 2 kV 100 kHz repetition rate	± 2 kV 100 kHz repetition rate
Surge voltages wire vs. wire	IEC 61000-4-5	± 0.5 kV, ± 1 kV	± 0.5 kV, ± 1 kV
Surge voltages wire vs. earth	IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2 kV	± 0.5 kV, ± 1 kV, ± 2 kV
Wire-conducted disturbances,	IEC 61000-4-6	3 V 0.15 MHz - 80 MHz 6 V in	3 V 0.15 MHz - 80 MHz 6 V in

Phenomenon	EMC basic standard	Immunity test level for professional health care facility	Immunity test level for areas of domestic health care
induced by high-frequency fields and connection to patients		ISM frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	ISM and amateur radio frequency bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Voltage drop	IEC 61000-4-11	0 % U_T ; 0.5 period at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % U_T ; 1 period and 70 % U_T ; 25/30 period monophasic: at 0 degrees	0 % U_T ; 0.5 period at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % U_T ; 1 period and 70 % U_T ; 25/30 period monophasic: at 0 degrees
Voltage interruptions	IEC 61000-4-11	0 % U_T ; 25/30 periods	0 % U_T ; 25/30 periods

**NOTICE**

U_T is the alternating voltage before application of the test level.

A.4.4 Test specifications for immunity and covers to HF wireless communication equipment

The device is intended for operation in an electromagnetic environment in which the HF disturbances are controlled. The customer or user of the device can thereby help to avoid electromagnetic interference by maintaining the minimum distance between portable and mobile HF telecommunication devices (transmitters) and the device, depending on the output power of the communication device.

Test frequency (MHz)	Frequency band (MHz)	Wireless service	Modulation	Max. power (W)	Distance (m)	Immunity test level (V/m)
385	380...390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27
450	430...470	GMRS 460 FRS 460	FM $\pm$ 5 kHz rise 1 kHz Sinus	2	0.3	28
710 745 780	704...787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
810 870 930	800...960	GSM 800/900 TETRA 800 iDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
1720 1845 1970	1700...1990	GSM 1800 CDMA 1900 GSM 1900	Pulse modulation 217 Hz	2	0.3	28

Test frequency (MHz)	Frequency band (MHz)	Wireless service	Modulation	Max. power (W)	Distance (m)	Immunity test level (V/m)
		DECT LTE Band 1, 3, 4, 25 UMTS				
2450	2400...2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240 5500 5785	5100...5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9

**NOTICE**

If necessary to achieve the immunity test levels, the distance between the transmitting antenna and the ME device or ME system can be reduced to 1 m. A testing distance of 1 m is allowed by IEC 61000-4-3.

A.5 Training

After having been instructed and having read the operating instructions, training on the device is not essential.

If there is a need, however, further training can be provided. Please contact PHYSIOMED ELEKTROMEDIZIN AG.


A.6 Settings

In the **Setting** menu, you can set the brightness of the screen and switch the dialog language.

You can access the basic settings menu from the welcome screen which is displayed as soon as you insert the therapy card.

How to Configure the Settings

- (1) Click the **SETTING**  symbol. The settings are displayed.

Settings

- (2) Click the selected language and rotate the **Data Selector <2>** to set the desired dialog language. Possible values are: **deutsch, english, español, français**
- (3) Click to confirm, and then click the brightness symbol.
- (4) Rotate the **Data Selector <2>** until the parameter has the desired value.
- (5) Click the  symbol to return to the menu.
The new settings are immediately active.

A.7 Technical Data

Treatment

Protection class according to IEC 60601-1	Internal power supply, Type BF
---	--------------------------------

Charging

Protection class according to IEC 60601-1	II
Input voltage	12 VDC
Input current	1.3 ADC

General Technical Data

CE label	according to regulation (EC) 2017/745 on medical devices
Class according to regulation (EC) 2017/745	Ila
Output voltage (max.)	400 Vs
Power supply	1 x 7.2V Li-Ion 24 Wh battery (fixed)
Output impedance	10 MOhm
Output frequency	5 ... 250 Hz
Duty cycle	1/4, 1/3, 1/2, 2/3, 3/4
Dimensions (W x H x D)	10.0 x 3.1 x 19.0 cm
Weight	0.5 kg
Degree of protection	IP 22: Protected against solid particles with a diameter ≥ 12.5 mm; protection against falling drip water

**Important**

Please note that a high storage temperature leads to increased self discharge of the battery!

Power Adapter

Type (to be used exclusively)	1 x 7,2V Li-Ion 24Wh (firmly installed)
Mains voltage	100 ... 240 VAC
Input current	0.35 A
Mains frequency	50 ... 60 Hz
Output voltage	7.2 VDC
Output current	1.3 ADC

A.7.1 Environment Conditions**Operation**

Ambient temperature	+5°C ... +40°C
Relative humidity	15% ... 90%
Air pressure	700 ... 1.060 hPa

Storage and Transport

Ambient temperature	-10°C ... +50°C
Relative humidity	10% ... 90%
Air pressure	500 ... 1,060 hPa

Appendix B Scope of Delivery and Accessories

Important

For safety reasons, the device is to be used exclusively with original accessories. Use of accessories by other manufacturers occurs at the risk of the user!

Please carefully read the instructions and warnings as well as the contraindications before first using the device.

B.1 Scope of Delivery

B.1.1 DEEP OSCILLATION® Personal Basic

DEEP OSCILLATION® Personal Basic is supplied with the following accessories:

Ref. no.	Designation	Quantity
00379	Applicator handhold	1
00261	Connection cable DEEP OSCILLATION®	2
00262	Connection cable grey for adhesive electrodes	1
00381	Oscillator head 5 cm	1
00589	PHYSIOPADS adhesive electrode for DEEP OSCILLATION®, set of 4	1
00348	Powder	1
55965	Power adapter Sinpro HPU15-105	1
00346	Special gloves size M	1
00372	Therapy card Basic	1
00382	Titanium neutral element	1
00395	Transportation bag	1
GA02040	Operating instructions (multilingual - DA, DE, EN, ES, FR, SV)	1

B.1.2 DEEP OSCILLATION® Personal Pro

DEEP OSCILLATION® Personal Pro is supplied with the following accessories:

Ref. no.	Designation	Quantity
00261	Connection cable DEEP OSCILLATION®	1
00262	Connection cable grey for adhesive electrodes	1
00371	Hand applicator set, consisting of: • 1x 00379 Applicator handhold	1

Ref. no.	Designation	Quantity
	2x 00381 Oscillator head, 5.0 cm 2x 00386 Oscillator head, 9.5 cm 1x 00261 Connection cable DEEP OSCILLATION®	
00589	PHYSIOPADS adhesive electrode DEEP OSCILLATION®, set of 4	1
00348	Powder	1
55965	Power adapter Sinpro HPU15-105	1
00346	Special gloves size M	1
00373	Therapy card Pro	1
00382	Titanium neutral element	1
00395	Transportation bag	1
GA02040	Operating instructions (multilingual - DA, DE, EN, ES, FR, SV)	1

B.2 Available Accessories

The following accessories are available for DEEP OSCILLATION® Personal:

Ref. no.	Designation
00379	Applicator handhold
00261	Connection cable DEEP OSCILLATION®
00262	Connection cable grey for adhesive electrodes
00393S	Membrane set Ø 1.5 cm (10 items)
00391S	Membrane set Ø 5 cm (10 items)
00392S	Membrane et Ø 9.5 cm (10 items)
00353	Monouse adhesive electrode (60 pcs.)
00387	Oscillator head 1.5 cm
00381	Oscillator head 5 cm
00386	Oscillator head, 9.5 cm
00396	Pin applicator
00348	Powder
55965	Power adapter Sinpro HPU15-105
00347	Special gloves size L
00346	Special gloves size M
00345	Special gloves size S
00420	Special gloves size XL
00355	Spiral lead black for adhesive electrodes
00356	Spiral lead universal
00338	Sweat band white (2 pcs.)
00382	Titanium neutral element

Ref. no.	Designation
00395	Transportation bag

B.3 CE Designation

Legal Note: The accessories are designated with the CE label (without number), where applicable.

The symbols have the following meaning:

Symbol	Meaning
	The accessory item is assigned to class 1 according to Annex VIII of the Council Regulation (EC) No. 2017/745 and complies with the corresponding requirements. Therefore it is labelled with the symbol on the left.
	Reference number
	Manufacturer
	Date of fabrication, year
	Retailer


Image similar

Appendix C Manufacturer's Recommendations

MANUFACTURER'S RECOMMENDATIONS
SAFETY REGULATIONS CONTROL
according to Medical Devices Regulation (MDR)

DEVICE: **DEEP OSCILLATION® Personal**
MANUFACTURER: PHYSIOMED ELEKTROMEDIZIN AG

The device has to undergo a safety regulation control every 12 months.

EXTENT:

- (1) Visual inspection of the device, accessories and accompanying papers
- (2) Function of controls and indicators
- (3) Functional testing of device and accessories
- (4) Curve shapes of output parameters
- (5) Output current at the ports
- (6) Electrical safety according to EN 62353, Substitute patient leakage current:
Limiting value according to EN 62353: 5 mA
Value first measured (NEW DEVICE): 0,05 mA

Appendix D Warranty

Please refer to the enclosed guarantee policy for terms of warranty.

Exceptions

Explicitly excluded from warranty are:

- Wear due to appropriate use or consumption of the product
- Accessory parts supplied with this product, which are worn or consumed due to appropriate use (e.g. batteries, electrodes, chipcards etc.)
- Products, which had been used, cleaned, stored or maintained inappropriately or against the terms of these operating instructions, as well as products that had been opened, repaired or modified by the purchaser or an authorized service partner
- Damage that occurred on the transportation route between the manufacturer and the customer, or between the service center and the customer, respectively
- Secondary damages caused by defects of this product (in this case demands due to product liability or other mandatory legal liability regulations might be claimed)

The original warranty time (beginning with the date of purchase) remains also in the case of repairs or after a full replacement and is not extended.

Index

A

accessories 30
 additional 31
 designation 32
 supplied 30
acute musculoskeletal pain 7
adverse events 8
ambient temperature 28
application 7

B

battery 13
brightness 27

C

card reader 11
cardiac arrhythmia 8
cardiac edema
 decompensated 8
cardiac insufficiency 8
CE label 28
charger socket 11
chronical obstructive pulmonary disease 7
claims 7
cleaning 23
compatibility
 electromagnetic 24
contraindications 8
control light 12
controls 10
 function 10
customer service 23

D

data selector 11
degree of protection 28
designation
 accessories 32
device description 7
device errors 18
device overview 4
dialog language 27
dimensions 28
disinfection 23
display 10
 brightness 27
disposal 23

E

economy mode 14
effects 7
EMC 24
environment conditions 29
error code 18
erysipelas 8
essential performance 6

F

fibromyalgia 7
function 7
function check 14

G

general notes 6

H

hamstring muscle tightness 7
hand applicator 15
hypersensitivity
 to electrostatic fields 8

I

implants
 electronic 8
indications 19
indicators 10
 function 10
individual programs 21
infections
 acute 8
inflammations
 acute 8
insufficiency
 vertebrobasilar 8
intended duration 9
intended frequency of contact 9
intended location 9
intended point of use 9
intended purpose 7
intended user 9
Intensity 11

L

lipedema 7
lymphedema 7

M

maintenance 23
malignant processes 8
manufacturer's recommendations 33
medical indications 7

membrane
 replace 15
modes of operation 7
modulation 28

O

output frequency 28
output impedance 28
output voltage 28

P

pacemakers 8
Parkinson's disease 7
pin applicator 16
ports 11
power adapter 13
power supply 11
precautions for use 9
pregnancy 8
protection class 28

R

rating plate 5
Raynaud phenomenon 7
risks 8

S

selftest 14, 18
service 23
service life 23
settings 27
shut-off 14
skin burns 7
skin diseases
 infectious 8
sleep mode 14
start-up 14

T

target patient group 7
technical data 28
therapeutic fields 7
therapy card 16
thromboses
 untreated 8
training 27
treatment 17
tuberculosis
 active 8

V

VBI

contraindication 8
venous diseases
 acute 8

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Copyright® PHYSIOMED ELEKTROMEDIZIN AG