LAVANID® - Wundspüllösung

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en

INSTRUCTIONS FOR USE

LAVANID®-wound irrigation solution

Read the instructions for use completely before using the product.

Use the product only for its intended purpose.

Please note the meanings of the symbols at the end of the instructions for use.

Note for patients:

This product is used by a doctor or other healthcare professional for the treatment of wounds. Patients

should be educated by a doctor or other healthcare professional before use. The type and duration of use is determined by the condition of the wound. If you notice any adverse effects, contact your doctor or health care professional.

DESCRIPTION AND INTENDED PURPOSE

This product is a sterile, ready-to-use wound irrigation solution for cleaning and moistening wounds. This product is suitable for acute wounds and for repeated and long-term wound treatment of chronic wounds. Chronic wounds are often covered with cell debris and dead (necrotic) tissue and may be colonised by germs. Wound coatings can lead to delayed wound healing and should be removed during wound treatment, for example by debridement. Gentle wound cleansing and wound moistening is necessary as an accompanying measure.

LAVANID®-wound irrigation solution consists of an isotonic salt solution, the so-called Ringer's solution, with 0.04 % polyhexanide for preservation. The preservative polyhexanide prevents the growth of microorganisms in wound dressings and in the solution. The chosen combination of Ringer's solution and polyhexanide has been well tolerated by the tissue in experimental and clinical studies. Wound treatment products preserved with polyhexanide are considered as products of choice for chronic wounds and burns. In a clinical study on chronic wounds, LAVANID®-wound irrigation solution was well tolerated. Yellowish wound coatings and mucous films containing microorganisms (biofilms) as well as unpleasant wound odours and signs of inflammation were reduced. The healing process is characterised by a decrease in wound size.

COMPOSITION

Sodium chloride, potassium chloride, calcium chloride, water for injection, polyhexanide.

FIELDS OF APPLICATION

LAVANID®-wound irrigation solution is applied:

- for local wound irrigation and wound cleansing
- for intraoperative wound irrigation during surgical procedures
- for diagnostic and therapeutic measures that require irrigation or moisturising
- for moistening cloths, tamponades and bandages

A post-treatment or rinsing off is not necessary.

CONTRAINDICATION

LAVANID®-wound irrigation solution must not be used:

- in the entire area of the central nervous system and the membranes of the brain and spinal cord (meninges)
- in the area of cartilage and joints. Risk of cartilage damage!
- in the middle ear or inner ear
- in the abdominal cavity (intraperitoneal)
- in case of known hypersensitivity (allergy) to one of the ingredients
- during pregnancy and breastfeeding

WARNINGS AND PRECAUTIONS

The product is sterile as long as the bottle is unopened and undamaged.

Do not use if the product shows signs of damage or if the packaging has been damaged or accidentally opened before use.

Do not resterilize or reprocess. This can lead to critical changes in performance, safety and handling of the product.

Do not use for injection or infusion.

Do not drink or swallow the product.

Only use clear solutions!

Do not mix with additives.

When handling and using sterile irrigation solutions, the principles of germ-free (aseptic) work must be observed.

INTERACTIONS WITH OTHER PRODUCTS

When used together with other wound treatment products, the instructions for use of these products should be checked for possible interactions with the ingredients of LAVANID®-wound irrigation solution. Due to incompatibilities with the ingredient polyhexanide, LAVANID®-wound irrigation solution must not be used together with anionic substances. These are, for example, substances that reduce the surface tension (anionic surfactants).

Direct contact of PVP iodine preparations with LAVANID®-wound irrigation solution in the same wound area should be avoided. Superficial skin disinfection with PVP-iodine preparations before a

surgical procedure and subsequent irrigation with LAVANID® - wound irrigation solution during the surgical procedure is possible.

ADVERSE EFFECTS

The ingredient polyhexanide can cause allergic reactions. In individual cases, acute reactions of the immune system (anaphylactic reactions) have been reported.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

REMOVAL OF THE SOLUTION

The stopper allows aseptic removal of the irrigation solution, e.g. with a spike. Alternatively, the screw cap can be removed to pour out the solution. Further instructions for removing the solution can be found at the end of these instructions for use.

STORAGE AND SHELF LIFE

The product must not be used after the expiry date. The expiry date is indicated on the label.

30°C
2°C Temperature limit.
Store between 2°C and 30°C.

Store in a dry place.

Do not use if the product is stored outside the recommended storage conditions.

Fragile, handle with care.

Keep out of reach of children.

Shelf life after opening:
2 months.

Due to the preservation, a shelf-life of two months is possible after the first opening of the bottle. A prerequisite for this is that the solution is protected from contamination. The bottle must be closed immediately after each use. The bottle opening must be protected from contamination. The bottle must not come into contact with the wound. Bottles

which have come into contact with the wound or have been contaminated in any other way must be disposed of.

LAVANID® - wound irrigation solution can be warmed to body temperature (approx. 37°C) immediately before use.

DISPOSAL INFORMATION

The disposal of unused products or waste material must be carried out in accordance with specific national regulations.

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SERAG-WIESSNER GmbH & Co. KG