

REF	Glucose	Eiweiß	PH-Wert	Urobilinogen	Blut	Nitrit	Keton	Bilirubin	Dichte	Leukozyten	Albumin
C1 12					x					x	
C1 11	x	x	x								
C1 13	x	x			x						
C1 15	x	x	x			x					
C1 42	x	x			x	x				x	
C1 19	x	x		x	x	x				x	
C1 29	x	x	x	x	x	x	x	x	x		x
C1 291	x	x	x	x	x	x	x	x	x		x
C1 1001	x	x	x	x	x	x	x	x	x	x	x
C1 1001M-TEST	x	x	x	x	x	x	x	x	x	x	x
C1 11544373	x	x	x	x	x	x	x	x	x	x	x
C1 1544039											x
C1 24	x										
C1 25	x					x					
C1 46						x					
C1 23	x	x	x		x	x	x			x	